



## Enhance Visualization and Detection

One look on your viewbox and you will see that the KODAK MIN-R EV Film System delivers outstanding advantages for diagnostic and screening mammography. Thanks to its unique film emulsion layer structure, it provides high contrast for better imaging of both the parenchyma and the periphery. Its slow-film/fast-screen combination yields fine grains, reducing noise and delivering sharp detail. It is a combination that improves visualization of faint, small objects throughout the breast—even in dense tissue. Supported by the Carestream commitment to women’s health through superior performance.

## ADVANTAGES AT A GLANCE

- Optimal image quality through advanced technology
- Fine grain for outstanding detail
- Bright, sharp results with “sparkle” and high contrast
- Developed by the world leader in screen-film mammography



# KODAK MIN-R EV Film System

## Optimal image quality with “sparkle” and fine detail

Unique emulsion grains produce a sharper toe, which yields whiter whites in glandular tissue. As a result, it helps reveal faint, small objects such as microcalcifications. Furthermore, it does so without losing contrast in fatty areas or breast periphery. So it improves image quality for all breast tissue types, giving you greater confidence in your interpretation.

## Screen-film characteristics<sup>1</sup>

Film	Screen	Processing Cycle	Relative Speed <sup>2</sup>	Contrast <sup>3</sup>		D-Max
				RP	EX II	
MIN-R EV	EV 150	Standard	150	4.5	4.7	>4.5
MIN-R EV	EV 190	Standard	190	4.5	4.7	>4.5


<sup>1</sup> Based on processing in KODAK X-OMAT Processors using seasoned KODAK X-OMAT EX II or KODAK RP X-OMAT Developer and KODAK RP X-OMAT LO Fixer at recommended temperatures.

<sup>2</sup> Relative speed determined from matched-density radiographs.

<sup>3</sup> Contrast measured as the average gradient between densities 0.25 and 2.00 above gross fog.

## Superb image results with a choice of screens

MIN-R EV Film works with the KODAK MIN-R EV 150 Screen for lower noise and the KODAK MIN-R EV 190 Screen for lower dose. The improved phosphor and coating structures of these screens deliver high MTF and brightness. Both screens are available mounted in a MIN-R 2 Cassette. The design of this cassette allows the film to be placed against the chest wall, enhancing chest-wall imaging.



## The integrated approach to superior performance

As an advocate for women’s health, particularly when it comes to the early detection of breast cancer, Carestream is your trusted partner. Our mammography films lead the industry, and our integrated products, service, and support deliver superior performance. Plus, you can count on Carestream expertise when it comes to quality control, troubleshooting, education, and more. We are dedicated to enabling you to provide diagnostic and screening mammography that’s reliable and affordable.

## Ordering information

### MIN-R EV Film

- 100-sheet packages, not interleaved
- 18 x 24 cm
- 24 x 30 cm

### MIN-R EV 150 Screen in MIN-R 2 Cassette

- 18 x 24 cm
- 24 x 30 cm

### MIN-R EV 190 Screen in MIN-R 2 Cassette

- 18 x 24 cm
- 24 x 30 cm

### KODAK Screen Cleaner

- 1 oz (473 ml) bottle

### MIN-R Screen Cleaner Wipes

- 1 container, 50 wipes

