



Required Fields

DIRECTVIEW software allows the Key Operator to choose from a default set of patient data fields and set the importance of any field. Each field can be set to act according to specific rules:

- **Normal Textbox** allows the technologist to accept information sent from the worklist or to enter data in the field if needed
- **Required to Submit** reminds the user to fill the field before they start imaging
- **Required to Deliver** only alerts the user that information is missing after an image has been accepted to be sent to an appropriately configured destination.
- **Read Only Textbox** will not allow the user to change information sent from the RIS.
- **Read Only Label** is a field that will display, but hides the information in it. For example, there might be information that is pertinent to the interpretation of the study but for privacy reasons doesn't display on the CR/DR console.

Additionally, the Key Operator can set the importance of any field.

For example, if the site wants the Tech ID to be entered for every study, the Tech ID field can be set as a field that is "Required to Submit".

This designation will remind the user to enter a Tech ID before they can start imaging.

This information could prove to be helpful when the need arises to determine who did a particular study.

For example, when doing a reject analysis or addressing infection precautions when an infectious patient is identified after they've been to Radiology

The pink layout buttons are the set of the fields the Key Operator has to choose from.


Layout Configuration Laurie Cesar  1:12:05 PM

Image Viewer Editors | Markers Editor | Tech Assist | Patient Names | Patient Worklist | Patient Input | Image Review | RPDES


Patient Last Name test	Patient First Name Ili	Patient Middle Name	Patient ID	Gender	Date Of Birth	Exam Information
Patient Location	Department Name	Referring Physician	Priority	Study Date 2/9/2015 2:45:21 PM		More Information
Procedure Name Scoliosis - 1 View	Procedure Code	Accession Number	Tech ID Laurie Cesar			More Image Data
View Name Leg Length - AP - Erect						

Required / Editable





- Normal Textbox
- Required To Submit

Hospital Default without Middle Name | Hospital Default with Middle Name | Hospital with Multi-Byte Name | Hospital Default without Middle Name with CR Integration

Hospital Default with Middle Name with CR Integration | Hospital with Multi-Byte Name with CR Integration

Close 

Layout

 Save Changes  Main Menu  Cancel 

In this example, a DR system is used to scan CR cassettes, so the layout needs to include a field for Cassette ID. It is missing from this layout.

Layout Configuration Laurie Cesar 23 1:12:19 PM

Image Viewer Editors | Markers Editor | Tech Assist | Patient Names | Patient Worklist | Patient Input | Image Review | RPDES

Patient Last Name test	Patient First Name Ili	Patient ID	Gender	Date Of Birth	Exam Information More Information More Image Data
Patient Location	Department Name	Referring Physician	Priority	Study Date 2/9/2015 2:45:21 PM	
Procedure Name Scoliosis - 1 View	Procedure Code	Accession Number	Tech ID Laurie Cesar		
View Name Leg Length - AP - Erect					

Required / Editable

- Normal Textbox
- Required To Submit
- Required To Deliver to Mandatory Destinations
- Read Only Textbox
- Read Only Label

Select New Layout

Save Changes | Main Menu | Cancel

The screen below shows that **Cassette ID** and **Orientation** fields are now available as a result of choosing the "Hospital Default without Middle Name **with CR Integration**" button.

The screenshot displays the 'Layout Configuration' window for user Laurie Cesar. The interface includes a top navigation bar with buttons for 'Image Viewer Editors', 'Markers Editor', 'Tech Assist', 'Patient Names', 'Patient Worklist', 'Patient Input', 'Image Review', and 'RPDES'. The main area contains a grid of input fields for patient and procedure data. The 'Cassette ID' field is highlighted with a yellow border. A 'Required / Editable' legend is open, showing options for 'Normal Textbox', 'Required To Submit', 'Required To Deliver to Mandatory Destinations', 'Read Only Textbox', and 'Read Only Label'. The bottom of the screen features a 'Home' icon, 'Save Changes', 'Main Menu', and 'Cancel' buttons.

Field	Value
Patient Last Name	test
Patient First Name	Ili
Patient ID	
Gender	
Date Of Birth	
Patient Location	
Department Name	
Referring Physician	
Priority	
Study Date	2/9/2015 2:45:21 PM
Procedure Name	Scoliosis - 1 View
Procedure Code	
Accession Number	
Tech ID	Laurie Cesar
Cassette ID	
View Name	Leg Length - AP - Erect
Orientation	

Required / Editable

- Normal Textbox
- Required To Submit
- Required To Deliver to Mandatory Destinations
- Read Only Textbox
- Read Only Label

Buttons: Home, Save Changes, Main Menu, Cancel, Select New Layout