



site planning

**Kodak**  
DirectView  
DR 9500 System

# Unprecedented DR system flexibility


The KODAK DIRECTVIEW DR 9500 System delivers dual-detector versatility with the simplicity and convenience of a single-detector system. A unique U-arm design moves the tube and detector around the patient in continuous alignment.

This provides remarkable freedom and mobility for patient positioning—and unprecedented flexibility for all horizontal, upright, and angled projections traditionally performed with dual-detector systems.

A single integrated operator console, ergonomic controls, and easy importation of HIS/RIS patient data further save examination time and effort. The same intuitive touch-screen interface found in all DIRECTVIEW DR and CR systems simplifies operation and reduces training requirements.

It all adds up to better workflow, higher productivity, and enhanced patient comfort—with the superb image quality and reliability you expect from the Kodak brand.

Figure A  
Minimum room size



# Three major subsystems

The uniquely flexible KODAK DIRECTVIEW DR 9500 System includes three major components:

- ▶ Ceiling-mounted U-arm with Bucky and x-ray tube assembly
- ▶ X-ray generator
- ▶ Integrated operator console

## Rail Length

Min: 6.56 ft (2 m) Max: 16.4 ft (5 m)

## Rail Length

Min: 6.56 ft (2 m) Max: 19.68 ft (6 m)




Figure B

The DR 9500 accommodates a variety of room sizes. Rails are available in 40 in. (1 m) increments.



The illustrations in this document serve only as a room size guide and to depict key subsystems of the DR 9500 System. Full requirements for the system are outlined in Kodak publication no. 6H3663, and should always be used for detailed planning purposes. The site specs include, among other information, complete room, electrical, phone, network, and environmental requirements. Detailed component dimensions and weights are also provided.

# Kodak DirectView

## DR 9500 System

### Specifications

#### Minimum size requirements

Room size:	Min. (L x W) 9.84 x 9.84 ft (3 x 3 m)
Ceiling height:	9.35 ft (2.85 m)
Control area wall height:	6.99 ft (2.13 m)
Door size:	(L x W) 84 x 34 in. (213 x 86 cm)
Hallway width:	60 in. (152 cm)
Rail length:	Min. (L x W) 6.56 x 6.56 ft (2 x 2 m) Max. (L x W) 19.68 x 16.4 ft (6 x 5 m)
Site power:	208/230/240V 20 A for Console and Positioner 400 VAC, 480 VAC, 3-phase for Generator

Key components	Width	Height	Depth	Weight
X-ray generator:	16 in. (40.6 cm)	48.7 in. (123.7 cm)	19.2 in. (48.8 cm)	265 lb (119 kg)
Wall chase:	20 in. (50.8 cm)	96 in. (243.8 cm)(adjustable)	7 in. (17.8 cm)	—
Operator console:	33.9 in. (86 cm)	38.2 in. (97 cm)	22.8 in. (58 cm)	185 lb (83 kg)

#### Digital image capture device

Detector type:	Cesium iodide
Active image:	17 x 17 in. (43 x 43 cm)
Detector element pitch:	143 micron
Dynamic range:	14-bit linear

#### X-ray generator

Output:	80 kW high-frequency output with digital feedback control
Range of output:	1000 mA at 80 kVp 800 mA at 100 kVp 630 mA at 125 kVp 500 mA at 150 kVp
Output:	65 kW high-frequency output with digital feedback control
Range of output:	800 mA at 80 kVp 630 mA at 100 kVp 500 mA at 125 kVp 400 mA at 150 kVp

#### Ceiling-mounted U-arm system

Variable SID:	40–79 in. (1–2 m)
---------------	-------------------

#### Multiplanar motions with electronic locks:

Bucky tilt:	+/- 45°
Tube rotation:	+ 210° to -20°
U-arm rotates:	+ 225° to -45° on horizontal axis
U-arm rotates:	+/- 270° on vertical axis
Vertical drive when horizontal:	15.75–98.42 in. (40–250 cm)

#### U-arm (external dimensions) for a SID of 79 in. (2 m):

Width:	101 in. (2.57 m)
Height:	41.5 in. (1.05 m)

#### Network capability

DICOM 3.0:  
Print service class  
Storage service class  
Multiple destination routing

Connectivity:  
KODAK DRYVIEW Laser  
Imagers and KODAK  
CARESTREAM PACS System  
third-party DICOM-  
compliant devices

Patient demographic and  
exam data:  
Through HIS/RIS with  
optional KODAK DIRECTVIEW  
DICOM Worklist Management Software

#### Regulatory

System complies with 21 CFR  
Chapter 1, Subchapter J Collimator:  
Class 2m laser product in  
accordance with IEC 60825-1



## Kodak

Licensed Product

Carestream Health, Inc.  
150 Verona Street  
Rochester, NY 14608

Carestream Health Canada Company  
6 Monogram Place  
Suite 200  
Toronto, Ontario M9R 0A1  
CANADA

#### More information

Learn more about the KODAK DIRECTVIEW DR 9500 System—and Carestream Health's flexible family of DR systems. Contact your Carestream Health representative, or call 1-877-865-6325, ext. 227.

→ [www.carestreamhealth.com](http://www.carestreamhealth.com)

Carestream   
HEALTH