

Carestream

PRICE LIST

2014

Film
Systems
Oncology
Laser Imaging
Radiology

Contents

Film Systems **3**

Mammography **7**

Oncology & Radiotherapy **8**

Laser Imaging Systems **9**

Service & Support **11**

Health & Safety & Environmental Issues **13**

Business Partners **14**

GREEN-SENSITIVE

CARESTREAM T-MAT G/RA Film

A high-contrast, green-sensitive film featuring outstanding resolution. Its pre-hardened emulsion allows it to be processed in CARESTREAM LE Plus Developer either in a standard processor or CARESTREAM X-OMAT RA Processor in the 38 second K/RA cycle. However, the film also retains the flexibility to be processed in Rapid or Standard cycles using CARESTREAM RP X-OMAT Chemicals. This film is the prime choice for general radiography and where high definition is needed such as extremity radiography. CARESTREAM T-MAT RA Films are designed for use with CARESTREAM LANEX Screens and X-OMAT Cassettes for maximum diagnostic accuracy with a minimum of radiation exposure.

Centimetres	Sheets/box	Cat. No.	Price (£)
TMG/RA not interleaved			
15 x 30	100	521 5603	161.72
18 x 24	100	835 3492	121.19
18 x 43	100	521 7740	210.01
24 x 30	100	151 0023	196.39
30 x 40	100	167 6469	320.50
35 x 35	100	830 1905	336.71
35 x 43	100	137 3125	406.69

CARESTREAM T-MAT L/RA Film

A low-contrast, green-sensitive film featuring outstanding resolution. Its prehardened emulsion allows it to be processed in CARESTREAM LE Plus Developer either in a standard processor or CARESTREAM X-OMAT RA Processor in the 38 second K/RA cycle. It is ideal for examinations where low contrast is needed, for example chest radiography and extremities where soft tissue imaging is required. It is the ideal film where exposure control is difficult. CARESTREAM T-MAT RA Films are designed for use with CARESTREAM LANEX Screens and X-OMAT Cassettes for maximum diagnostic accuracy with a minimum of radiation exposure.

Centimetres	Sheets/box	Cat. No.	Price (£)
TML/RA not interleaved			
15 x 30	100	521 2329	161.72
18 x 24	100	841 0789	121.19
18 x 43	100	831 0286	210.01
24 x 30	100	170 9039	196.39
30 x 40	100	155 8469	320.50
35 x 35	100	121 9849	336.71
35 x 43	100	181 4631	406.69

BLUE-SENSITIVE

CARESTREAM MXBE Film

MXBE film is a blue sensitive, general purpose film for all applications and can be used with all blue and ultraviolet screens.

Centimetres	Sheets/box	Cat. No.	Price (£)
not interleaved			
18 x 24	100	771 0783	128.23
18 x 43	100	771 1476	221.74
30 x 40	100	771 0809	335.54
35 x 35	100	771 0817	355.96
235 x 45	100	771 0825	430.05

CARESTREAM EKTASCAN B/RA Film

Can be processed in CARESTREAM X-OMAT LE Plus Chemicals in the RA cycle. If converting from EKTASCAN NB Film, multifformat camera settings will need to be adjusted.

Centimetres (inches)	Sheets/box	Cat. No.	Price (£)
20.3 x 25.4 (8 x 10)	100	157 7501	245.87
35 x 43 (14 x 17)	100	817 4385	735.03

CARESTREAM X-OMAT Duplicating RA Film

A duplicating film with a straight line contrast to ensure excellent copies. It can be processed across all processing cycles. It may be handled under a suitable safelamp fitted with a CARESTREAM 6B or GBX-2 Safelight Filter.

Centimetres (inches)	Sheets/box	Cat. No.	Price (£)
18 x 24	100	114 1159	212.27
24 x 30	100	149 4418	353.96
30 x 40	100	172 2388	578.70
35 x 35	100	195 8156	607.72
35 x 43	100	837 3805	735.03
20.3 x 25.4 (8 x 10)	100	811 9117	258.92

not interleaved, packed as 2 x 50 sheets per box

CARESTREAM X-OMAT LE Plus Developer (for rapid and normal processing cycles)

Recommended for use primarily in modern rapid cycle processors, this chemical contains no glutaraldehyde. It requires no starter solution and comes as a single-part concentrate which facilitates fast and accurate mixing. It features enhanced batch-to-batch consistency and produces less odour than conventional chemicals and is, therefore, environmentally more acceptable. It has a lower oxidation rate which permits reduced replenishment rates and a longer useful life, all of which results in increased economy.

Description	Cat. No.	Price (£)
2 x 20 litres	526 9055	125.80

CARESTREAM RP X-OMAT EXII Developer

This new developer is non-corrosive, lower in odour, BOD and COD and the extremely low glutaraldehyde content is now complexed to make it more environmentally friendly.

Description	Cat. No.	Price (£)
To make 40 litres (2 x 20 litres)	527 4394	114.25

CARESTREAM RP X-OMAT Developer Starter

Description	Cat. No.	Price (£)
3 litres (6 x 500ml)	523 9389	46.40

CARESTREAM X-OMAT LE Plus Fixer (rapid and normal processing cycles)

Recommended for use primarily in modern rapid cycle processors, this single part fixer contains no hardener and can be used in processors with in-line mixing.

Description	Cat. No.	Price (£)
2 x 20 litres	528 4484	70.57

CARESTREAM RP X-OMAT LO-4 Fixer and Replenisher

This product is designed for automatic processing using a cycle of between 45 seconds and 3 minutes 30 seconds. It may be mixed manually, or using a CARESTREAM Automixer.

Description	Cat. No.	Price (£)
To make 40 litres (2 x 20 litres)	522 4381	60.06

CARESTREAM READYMATIC Developer and Replenisher

A ready-to-use solution intended for processing X-ray films in a CARESTREAM X-OMAT Processor, and those designed for dental X-ray film.

Description	Cat. No.	Price (£)
10 litres (2 x 5 litres)	502 3866	48.46

CARESTREAM READYMATIC Fixer and Replenisher

A ready-to-use solution for use with CARESTREAM READYMATIC Developer and Replenisher.

Description	Cat. No.	Price (£)
10 litres (2 x 5 litres)	502 3874	38.37

CARESTREAM X-OMAT MINI KITS

The CARESTREAM X-OMAT Mini Kit is formulated to optimise image quality, enhance your working environment and comply with environmental regulations. The CARESTREAM X-OMAT Developer and Replenisher, when used with CARESTREAM X-OMAT Fixer and Replenisher provides high quality images and consistent performance. The CARESTREAM X-OMAT Mini Kit is recommended for all CARESTREAM diagnostic films and is compatible with other manufacturers' films.

Developer: A convenient two-part liquid and developer concentrate provides 2.5 litres of developer solution

Fixer: The single part fixer concentrate provides 2.5 litres of fixer and replenisher.

Description	Cat. No.	Price (£)
CARESTREAM X-OMAT Developer & Replenisher Mini Kit	521 0026	102.54
CARESTREAM X-OMAT Developer & Replenisher Kit	521 0034	55.26

Cassettes & Screens

FOR MEDICAL IMAGING

CARESTREAM X-OMAT Cassette C-1 with LANEX Screens

A durable lightweight cassette, supplied complete with intensifying screens. The highly advanced materials used ensure the lowest possible dose. Its curved sides expel air upon closing and provide precise film/screen contact. All sizes have a single latch for easy opening and a built-in plastic moulded light-lock to ensure light-proofing.

Centimetres	Cat. No.	Price (£)
Supplied with a pair of CARESTREAM LANEX Fine Screens		
18 x 24	181 4086	438.02
18 x 43	145 8967	581.99
24 x 30	885 0398	537.12
30 x 40	173 8657	766.13
35 x 35	899 9948	766.13
35 x 43	180 9854	900.86
Supplied with a pair of CARESTREAM LANEX Medium Screens		
18 x 24	828 9647	438.02
18 x 43	144 5527	581.99
24 x 30	877 1925	537.12
30 x 40	867 8377	766.13
35 x 35	859 7320	766.13
35 x 43	196 9898	900.86
Supplied with a pair of CARESTREAM LANEX Regular Screens		
18 x 24	129 0253	438.02
18 x 43	896 4256	581.99
24 x 30	859 4384	537.12
30 x 40	186 6433	766.13
35 x 35	828 9506	766.13
35 x 43	181 5018	900.86
15 x 30 (windowless)	850 4706	on request
Supplied with a pair of CARESTREAM LANEX Fast Screens		
18 x 24	898 0138	438.02
18 x 43	101 1642	581.99
24 x 30	836 0869	537.12
30 x 40	154 7017	766.13
35 x 35	827 8491	766.13
35 x 43	129 0428	900.86

PLEASE NOTE: The above CARESTREAM X-OMAT Cassettes represent those varieties and sizes ordered most frequently. Other varieties are available, as are replacement CARESTREAM X-OMATIC Screens. Please enquire through your local Business Partner if the items you require are not listed.

CARESTREAM Screen Cleaner

A special cleaner intended for use with all CARESTREAM X-OMATIC, LANEX, EC and MIN-R screens. The CARESTREAM screen cleaner is specially formulated to optimise the superb image quality of the CARESTREAM Screen Systems and helps eliminate dust and dirt that can cause minus density artifacts on the film.

Description	Cat. No.	Price (£)
16 US fl oz (473ml)	845 4977	40.01

MAMMOGRAPHY FILM SYSTEMS

CARESTREAM MIN-R EV Film System

A higher contrast, dual emulsion film, compatible with either CARESTREAM MIN-R EV 150 screens or Carestream MIN-R EV 190 screens. Improved resolution and compatibility with glutaraldehyde-free developers.

MIN-R EV Film (Centimetres)	Sheets/boxes	Cat. No.	Price (£)
18 x 24	100	892 5356	216.53
24 x 30	100	185 1609	368.03
MIN-R 2 Cassette supplied with a single CARESTREAM			
MIN-R EV 150 Screen			
18 x 24		111 7035	418.22
24 x 30		828 3061	501.86
MIN-R 2 Cassette supplied with a single CARESTREAM			
MIN-R EV 190 Screen			
18 x 24		138 7414	483.29
24 x 30		880 6325	576.21

CARESTREAM MIN-R 2000 PLUS Film System

This film is compatible with either CARESTREAM MIN-R 2000 PLUS or CARESTREAM MIN-R 2190 Screens, and is designed specifically for standard processing.

MIN-R 2000 PLUS Film (Centimetres)	Sheets/box	Cat. No.	Price (£)
18 x 24	100	604 4655	206.33
24 x 30	100	604 4663	351.28
MIN-R 2 Cassette supplied with a single CARESTREAM			
MIN-R 2000 Screen			
18 x 24		810 4101	408.92
24 x 30		899 9492	445.11
MIN-R 2 Cassette supplied with a single CARESTREAM			
MIN-R 2190 Screen			
18 x 24		126 0579	468.40
24 x 30		891 0853	563.18

CARESTREAM MIN-R Screen Wipes

The New CARESTREAM MIN-R Screen Wipes are designed to ensure optimum screen cleanliness and minimise screen generated artifacts for best image quality.

Description	Cat. No.	Price (£)
50 wipes per cylinder		
Wipes 10 x10 cm lint free, pre-moistened	528 0094	19.57

IMAGING STANDARD FOR PORTAL IMAGING

Please note: For more information, or to place an order please contact:

Oncology Imaging Systems Ltd.

Tel: 01825 840633

Fax: 01825 841496

Email: sales@oncologyimaging.com

Laser Imaging Systems

FILMS AND PAPERS FOR LASER PRINTERS

CARESTREAM DRYVIEW CHROMA Imager

Delivers diagnostic quality images on film and referral prints on paper in a choice of high-quality greyscale or brilliant full colour.

CARESTREAM DRYVIEW CHROMA Medical Film and Papers for CARESTREAM DRYVIEW CHROMA Printer

CARESTREAM DRYVIEW CHROMA Medical Film

Centimetres (inches)	Sheets/box	Cat. No.	Price (£)
20 x 25	500 (5 x 100 shts)	604 4580	163.96
20 x 25	10 x 10 starter packs	604 4614	on request
25 x 30	500 (5 x 100 shts)	604 4572	242.26
28 x 35	500 (5 x 100 shts)	604 4564	301.45
35 x 43	500 (5 x 100 shts)	604 x 4556	448.62

CARESTREAM DRYVIEW CHROMA Glossy Paper

Centimetres (inches)	Sheets/box	Cat. No.	Price (£)
28 x 75m	75m roll	604 4937	on request
29.7 x 75m	75m roll	604 4929	177.33
21.5 x 28	10 x 10 starter packs	604 4606	on request
21 x 29.7	10 x 10 starter packs	604 4598	on request
21.5 x 28	500 (5 x 100 shts)	604 4549	on request
21 x 29.7	500 (5 x 100 shts)	604 4531	152.89
29.7 x 42	500 (5 x 100 shts)	604 4523	282.88

CARESTREAM DRYVIEW CHROMA Mat Paper

Centimetres (inches)	Sheets/box	Cat. No.	Price (£)
28 x 100m	100m roll	604 5025	on request
29.7 x 100m	100m roll	604 4945	157.63
21.5 x 28	750 (5 x 150 shts)	604 4515	on request
21 x 29.7	750 (5 x 150 shts)	604 4507	164.02
29.7 x 42	700 (5 x 150 shts)	604 4499	244.26

CARESTREAM DRYVIEW Ink Cartridges

Description	Cat. No.	Price (£)
Photo Black Ink	655 7813	92.87
Cyan Ink	655 7821	92.87
V Magenta Ink	655 7839	92.87
Yellow Ink	655 7847	92.87
Light Cyan Ink	655 7854	112.87
VL Magenta Ink	655 7862	112.87
Light Black Ink	655 7870	112.87
Mat Black Ink	655 7888	112.87
L Light Black Ink	655 7896	92.87
Orange Ink	655 7904	92.87
Greene Ink	655 7912	92.87

Laser Imaging Systems

FILMS AND PAPERS FOR LASER PRINTERS

CARESTREAM DRYVIEW PLUS Laser Films

A high definition dry imaging film suitable for use in CARESTREAM DRYVIEW Laser Imaging Systems.

Centimetres (inches)	Sheets/box	Cat. No.	Price (£)
For DRYVIEW 8100/8150/6800			
35 x 43 (14 x 17)	*500 (4 x 125 shts)	106 9590	3,065.52
For DRYVIEW 8900/6800			
35 x 35 (14 x 14)	*500 (4 x 125 shts)	190 9746	2,504.56
For DRYVIEW 8200/8900/6800			
27.9 x 35.6 (11 x 14)	*500 (4 x 125 shts)	149 5126	2,030.24
For DRYVIEW 8500			
27.9 x 35.6 (11 x 14)	*500 (4 x 125 shts)	163 1944	2,030.24
For DRYVIEW 8900/6800			
25.5 x 30.6 (10 x 12)	*500 (4 x 125 shts)	159 2229	1,569.11
For DRYVIEW 8900/6800			
20.3 x 25.4 (8 x 10)	*500 (4 x 125 shts)	149 1588	1,079.15

(*Minimum order 4 x 125 sheets - invoiced per 125 sheets †Minimum order 5 x 100 sheets - invoiced per 100 sheets)

For CARESTREAM DRYVIEW 5800

Description	Sheets/box	Cat. No.	Price (£)
For DRYVIEW DVB + CE			
20.3 x 25.4 (8 x 10)	500 (5 x 100 shts)	102 7697	1,079.15
For DRYVIEW DVB + CE			
25.5 x 30.6 (10 x 12)	500 (5 x 100 shts)	841 4666	1,569.11
For DRYVIEW DVB + CE			
27.9 x 35.6 (11 x 14)	500 (5 x 100 shts)	863 5294	2,030.24
For DRYVIEW DVB + CE			
35 x 35 (14 x 14)	500 (5 x 100 shts)	843 3518	2,504.56
For DRYVIEW DVB + CE			
35 x 43 (14 x 17)	500 (5 x 100 shts)	119 6740	3,065.52

CARESTREAM DRYVIEW DVM Film

CARESTREAM DRYVIEW DVM Film with improved contrast and image tone is optimised for mammography applications and is the only dry laser film available for digital mammography. It has high maximum density, so subtle shading remains crystal clear with lifetime archivability. Daylight packaging provides greater ease of use.

For CARESTREAM DRYVIEW 8610

Description	Sheets/box	Cat. No.	Price (£)
20.3 x 25.4 (8 x 10)	100	864 0823	226.64

For CARESTREAM DRYVIEW 8900/6800

Description	Sheets/box	Cat. No.	Price (£)
20.3 x 25.4 (8 x 10)	500 (4 x 125 shts)	145 6516	1,133.21
25.5 x 30.6 (10 x 12)	500 (4 x 125 shts)	808 4485	1,644.73

For CARESTREAM DRYVIEW DVM PLUS for mammo imaging (high Dmax = 4.0)

Description	Sheets/box	Cat. No.	Price (£)
DryView DVM+			
20.3 x 25.4 (8 x 10)	500 (5 x 100 shts)	141 1156	1,133.21
DryView DVM+			
25.5 x 30.6 (10 x 12)	500 (5 x 100 shts)	828 5090	1,647.03

CARESTREAM HEALTH EQUIPMENT MAINTENANCE AGREEMENTS

Equipment down time is expensive, causing problems for the care of patients and the morale of your team. A Carestream Health Equipment Maintenance Agreement offers you peace of mind in the knowledge that both your business requirements and operational needs are met.

All levels of Carestream Health Equipment Maintenance Agreements share the following benefits:

- Unlimited emergency support within normal hours and parts cover
- Access to technical assistance over the telephone
- Preventative maintenance to keep equipment at peak performance
- 'Best endeavour' cover at weekends
- One convenient annual invoice to cut administration

Extend your confidence

Carestream Health's flexible approach to service allows each Equipment Maintenance Agreement to be individually designed and built around your specific service needs and expectations. Choose from a multitude of offers and from a wide range of options to find the best package for you and be confident that we will be there to support you in your aims and objectives.

Tailor support from a range of choices including:

- 4, 8 or 16-hour response time*
- Extended working day 9.00am-9.00pm*
- Saturday mornings
- On-site or remote diagnostics
- Different levels of support for critical and non-critical equipment
- Longer term packages
- 24/7 cover for total confidence*

*Dependent on postal district—price on application.

In addition, choose from a range of additional services to really customise your package and add particular value to your operation.

- Detailed consultation prior to new investments
- Project design and management
- Workflow analysis
- DICOM V & V equipment compatibility and assessments
- Proactive remote diagnostics
- PACS/RIS integration and support
- Key operator training
- Network assessment, design and implementation
- Wireless detector drop waiver
- Ultrasound probe repair/replacement service

Whatever your needs, we can design a service package to suit you. Talk to your local Carestream Health Account Manager for a full costing based on meeting your individual requirements.

For easy reference keep details of your account manager and telephone number here:

Contact numbers

Equipment Maintenance Administration

0870 600 0245

Service and Support

UK, Ireland and Channel Islands

0870 600 0245

BEHIND THE SCENES

Whenever you need help, Carestream Health's extensive service infrastructure means that we can call on local, national, European and global expertise as appropriate. This depth of resources truly sets Carestream Health Service and Support apart from its competitors.

At your side when you need us

Our base of over sixty highly trained engineers can work confidently across both conventional and digital environments. Investment in training of engineers has always been a top priority and Carestream Health utilises independent training programmes to achieve the highest level of training across the industry.

When problems arise our dedicated Call Centre team moves quickly to assess the problem and redirect your call accordingly.

If an engineer visit is required, their detailed knowledge will allow the appropriate parts to be picked in advance cutting out a repeat visit.

Specialist help over the telephone

The Technical Assistance Centre (TAC) can give you immediate technical advice from specialists over the telephone. When an equipment fault is raised, the TAC engineer can draw information from our fault database and suggest ways of solving the problem remotely and if necessary, order emergency spare parts for our local engineer to fit. In some cases, dial-in diagnostics may be possible and TAC can remotely repair the equipment. When on-site, they can also work with remote support from the TAC for additional technical assistance. Either way TAC can solve problems faster and improve your up-time.

HEALTH & SAFETY & ENVIRONMENTAL ISSUES

The labels and instructions provided with CARESTREAM Products normally give general working and health, safety and environmental information. Where necessary, they will specifically mention any constituents of health and safety and environmental concern and provide precautionary information relating to the safe handling, use and disposal of the product. Safety Data Sheets for CARESTREAM Chemical Products and any additional information you or your staff may need are available from:

your local Business Partner

or visit the Carestream Health website on

www.carestream.co.uk

or write to:

**Carestream Health UK Ltd.
Product Safety Advisory Service
One Park Lane
Hemel Hempstead
Herts. HP2 4YJ**

**Carestream Health UK Ltd.
One Park Lane
Hemel Hempstead
Herts. HP2 4YJ**

**Tel: 0870 600 0245
Fax: 01442 263055
Website: www.carestream.co.uk**

Carestream Health UK Ltd.
Registration No. 5942397 London
Registered office:
Carestream Health UK Ltd.
One Park Lane
Hemel Hempstead
Herts. HP2 4YJ

Business Partners

Carestream Health supplies its products through two Business Partners who cover all areas of the United Kingdom and Ireland. These specialist companies are dedicated to providing Carestream Health customers with the highest quality service at the best value for money.

Please place all orders for consumables through your appointed Business Partner (see locations on the map below) and any equipment and service requirements via Carestream Health's Business Support Team.

For easy reference please keep your account number here:

Carestream Health's Business Support Team can be contacted on

0870 600 0245



Carestream Health UK Ltd.

One Park Lane
Hemel Hempstead
Herts. HP2 4YJ

Tel: 0870 600 0245
Fax: 01422 263055

Registered office:

Carestream Health UK Ltd.
One Park Lane
Hemel Hempstead
Herts. HP2 4YJ
Registration No. 5942397 London

Carestream

www.carestream.co.uk

CARESTREAM and all product names displayed in capitals are trademarks of Carestream Health.

Because of our constant endeavour to improve quality and design, modifications may be made from time to time.
Any product, particularly one of an experimental nature or available only to special order, may be discontinued without notice.
Details of stock availability and prices given in this publication are subject to change without notice.
All weights and dimensions are approximate.