



Security Bulletin - Digital Capture-Oncology Security Bulletin

As of 30-December-2008

Contact Carestream Health, Inc. Service to make arrangements for your equipment

Color Key Code:

Unacceptable risk

Low risk, fix will be in next major release or rollup patch

N/A or N/A due to hardening, but may be applicable to other products

Note: any patches not on this report are deemed not applicable to Carestream products.

Product Susceptibility Analysis			KROS		
Vulnerability ID	Vulnerability Description	Date Analyzed	5.1	Dosimetry 1.0 / 1.1	6.0 / 6.1 / 6.2
MS08-078	Security Update for Internet Explorer (960714)	12/18/2008	High Risk	High Risk	High Risk
MS08-073	Cumulative Security Update for Internet Explorer (958215)	12/12/2008	High Risk	High Risk	High Risk
MS08-070	Vulnerabilities in Visual Basic 6.0 Runtime Extended Files (ActiveX Controls) Could Allow Remote Code Execution (932349)	12/9/2008	N/A	N/A	N/A
MS08-067	Vulnerability in Server Service Could Allow Remote Code Execution (958644)	10/24/2008	High Risk	High Risk	High Risk
MSSA956391	Cumulative Security Update of ActiveX Kill Bits	10/15/2008	Low Risk	Low Risk	Low Risk
MS08-066	Vulnerability in the Microsoft Ancillary Function Driver Could Allow Elevation of Privilege (956803)	10/15/2008	N/A	Low Risk	Low Risk
MS08-064	Vulnerability in Virtual Address Descriptor Manipulation Could Allow Elevation of Privilege (956841)	10/15/2008	Low Risk	Low Risk	Low Risk
MS08-063	Vulnerability in SMB Could Allow Remote Code Execution (957095)	10/15/2008	Low Risk	Low Risk	Low Risk

Product Susceptibility Analysis			KROS		
Vulnerability ID	Vulnerability Description	Date Analyzed	5.1	Dosimetry 1.0 / 1.1	6.0 / 6.1 / 6.2
MS08-062	Vulnerability in Windows Internet Printing Service Could Allow Remote Code Execution (953155)	10/15/2008	Low Risk	Low Risk	Low Risk
MS08-061	Vulnerabilities in Windows Kernel Could Allow Elevation of Privilege (954211)	10/15/2008	Low Risk	Low Risk	Low Risk
MS08-058	Cumulative Security Update for Internet Explorer (956390)	10/15/2008	Low Risk	Low Risk	Low Risk
MSSA953839	Cumulative Security Update of ActiveX Kill Bits	8/13/2008	Low Risk	Low Risk	Low Risk
MS08-050	Vulnerability in Windows Messenger Could Allow Information Disclosure (955702)	8/13/2008	N/A	N/A	N/A
MS08-049	Vulnerabilities in Event System Could Allow Remote Code Execution (950974)	8/13/2008	Low Risk	Low Risk	Low Risk
MS08-047	Vulnerability in IPsec Policy Processing Could Allow Information Disclosure (953733)	8/13/2008	N/A	N/A	N/A
MS08-045	Cumulative Security Update for Internet Explorer (953838)	8/13/2008	Low Risk	Low Risk	Low Risk
MS08-037	Vulnerabilities in DNS Could Allow Spoofing (953230)	7/10/2008	Low Risk	Low Risk	Low Risk