



Take Mobile Imaging to a New Level— with X-Factor Thinking.

With the Universal Mobile retrofit kit, you can upgrade portable analog X-ray systems from almost any manufacturer to the power of DR wireless technology. Easy to install, this innovative kit helps you maximize your investment in your current mobile units.

The Universal Mobile Retrofit Kit is a non-invasive, flexible retrofit package. Self-powered and rechargeable, the DRX detector transfers images instantaneously, so they can be quickly transferred to physicians for diagnosis and treatment planning. And, since you won't have to contend with tethers or cables during examinations, this is a system you'll use with ease and convenience in any healthcare facility.

ADVANTAGES AT A GLANCE

- Quickly and cost-effectively convert existing analog mobile systems to DR.
- Easy to use, minimal training necessary.
- Wireless capability improves productivity and workflow.
- Self-powered and rechargeable
- Available with either DRX-1 or more dose sensitive DRX-1C detector.



CARESTREAM DRX-Transportable System/Universal Mobile

It fits, it's fast, it's cost effective.

- Upgrade existing equipment to DR in a matter of hours
- Easy to learn, easy to use
- Wireless capability lets you go anywhere
- Multiple tablet PC mounting options to fit your mobile X-ray system
- Electronics fit in the existing cassette bin



Product Specifications

Key Components:

- Thin electronics box 454 x 368 x 48 mm (17.9 x 14.5 x 1.9 in.) containing the wireless access point, ethernet switch, generator interface board, and power supply (fits into the existing cassette bin)
- Electronics box uses the same type of battery as the DRX detector
 - Battery charges while plugged in
- Motion Computing™ tablet PC
 - I7 Processor / 12.1 in. touch-screen / 1.6 kg (3.6 lbs)
 - Rugged design that meets MIL-STD-810G and IP52
 - Dual Hot Swappable battery
- Mounting hardware to mount tablet PC to system



X-Factor Thinking

The X-Factor makes the most of your capital spending with configurable, modular components powered by the portable DRX detector. Use the detector on your retrofitted mobile units on your morning rounds... and later, take it to a busy DRX-Evolution room for peak patient volumes. When the X-Factor lets you scale your solution to fit your workflow and budget, the possibilities are endless.

Trusted Partner – CARESTREAM Smart Services

At Carestream, we'll be there when you need us. Our medical imaging solutions are backed by our worldwide support team of product and industry specialists—experts who work closely with you to keep your equipment up and running at peak performance. By delivering a comprehensive scope of service offerings, we provide complete peace of mind over the life of your equipment, so you can focus on what's most important—your patients.



www.carestream.com

