

## Kodak RP X-OMAT LO FIXER AND REPLENISHER

### Working Solution (diluted solution)

Contains:	HMIS CODE
Water (007732-18-5)	Health - 1
Ammonium thiosulfate (007783-18-8)	Flammability - 1
Sodium thiosulfate (007772-98-7)	Reactivity - 0
	Personal Protection-A

Low hazard for recommended handling. In case of contact, wash with plenty of water. For additional information, see Material Safety Data Sheet (MSDS).



Concentrates (not diluted solution) made by  
Eastman Kodak Company  
Rochester, New York 14650

Product Information: 800-328-2910  
Emergency Information: 585-722-5151  
TM: Kodak and X-OMAT

## Kodak RP X-OMAT DEVELOPER REPLENISHER

### Working Solution (diluted solution)

Contains:	HMIS CODE
Water (007732-18-5)	Health - 2
Potassium Sulfite (010117-38-1)	Flammability - 0
Hydroquinone (000123-31-9)	Reactivity - 0
Potassium acetate (000127-08-2)	Personal Protection - C
Glutaraldehyde bis (potassium bisulfite) (068310-08-7)	

WARNING! MAY BE HARMFUL IF SWALLOWED. CAUSES EYE IRRITATION. MAY CAUSE ALLERGIC SKIN REACTION. IN CASE OF CONTACT, WASH WITH PLENTY OF WATER.

For additional information, see Material Safety Data Sheet (MSDS)



Concentrate (not diluted solution) made by  
Eastman Kodak Company  
Rochester, New York 14650

Product Information: 800-326-2910  
Emergency Information: 585-722-5151  
TM: Kodak and X-OMAT

## Kodak RP X-OMAT LO FIXER AND REPLENISHER

### Working Solution (diluted solution)

Contains:	HMIS CODE
Water (007732-18-5)	Health - 1
Ammonium thiosulfate (007783-18-8)	Flammability - 1
Sodium thiosulfate (007772-98-7)	Reactivity - 0
	Personal Protection-A

Low hazard for recommended handling. In case of contact, wash with plenty of water. For additional information, see Material Safety Data Sheet (MSDS).



Concentrates (not diluted solution) made by  
Eastman Kodak Company  
Rochester, New York 14650

Product Information: 800-328-2910  
Emergency Information: 585-722-5151  
TM: Kodak and X-OMAT

## Kodak RP X-OMAT DEVELOPER REPLENISHER

### Working Solution (diluted solution)

Contains:	HMIS CODE
Water (007732-18-5)	Health - 2
Potassium Sulfite (010117-38-1)	Flammability - 0
Hydroquinone (000123-31-9)	Reactivity - 0
Potassium acetate (000127-08-2)	Personal Protection - C
Glutaraldehyde bis (potassium bisulfite) (068310-08-7)	

WARNING! MAY BE HARMFUL IF SWALLOWED. CAUSES EYE IRRITATION. MAY CAUSE ALLERGIC SKIN REACTION. IN CASE OF CONTACT, WASH WITH PLENTY OF WATER.

For additional information, see Material Safety Data Sheet (MSDS)



Concentrate (not diluted solution) made by  
Eastman Kodak Company  
Rochester, New York 14650

Product Information: 800-326-2910  
Emergency Information: 585-722-5151  
TM: Kodak and X-OMAT

### ATTACH LABELS TO REPLENISHER AND/OR PROCESS MACHINE TANKS

**DANGER!**

**DO NOT ADD CLEANING AGENTS TO TANKS  
UNLESS RECIRCULATION LINES AND TANKS HAVE  
BEEN DRAINED AND RINSED WITH WATER**

**Mixing Cleaning Agents with Photochemicals May  
Liberate Hazardous gases.**

For processor cleaning methods, call 1-800-328-2910

Read the Material Safety Data Sheet for chemicals used.

*These labels are provided to assist you in complying  
with Federal OSHA Hazard Communication Standard  
29 CFR 1910.1200*