


MORE INFORMATION. HIGHER CONFIDENCE. BETTER CARE.

Visualize substantially more diagnostic information for
Neonatal and Pediatric Applications




With CARESTREAM INSIGHT Pediatric Imaging Systems, the radiologist can visualize substantially more diagnostic information than is available with conventional imaging systems, due to the extremely low noise and high resolution. At the same time, improved diagnostic information and excellent exposure tolerance mean fewer repeat examinations. Three systems offer a choice of speed, yielding various levels of detail for a variety of diagnostic requirements.

DIAGNOSTIC INFORMATION FOR ALL PEDIATRIC:

CARESTREAM INSIGHT Pediatric Imaging

- CARESTREAM INSIGHT Pediatric Film
- CARESTREAM INSIGHT Pediatric Screens
- Speed: 300


This system delivers the high spatial resolution needed to differentiate neonatal anatomical structures and catheters – without the increased noise that can obscure certain disease patterns. Low noise is especially valuable in the diagnosis of hyaline membrane disease in premature newborns.



CARESTREAM INSIGHT Pediatric Detail Imaging

- CARESTREAM INSIGHT Pediatric Film
- KODAK InSight Pediatric Detail Screens
- Speed: 230

This system discloses subtle breakdowns in the trabecular pattern of the bone which occurs when the normal cellular appearance becomes disrupted by metabolic changes that precede cancer formation. As a result, potentially serious bone conditions, such as cancer, may be detected and treated earlier.



EXAMINATIONS APPLICATIONS:

CARESTREAM INSIGHT Pediatric Ultra-Detail Imaging

- CARESTREAM INSIGHT Pediatric Film
- KODAK InSight Pediatric Ultra-Detail Screens
- Speed: 80

Ultra-detail imaging is essential for differential diagnosis of trauma vs. traumatic abuse. Periosteal reaction, metaphyseal fractures, callus formation, and other cortical and trabecular variations may provide the only evidence of non-accidental injury.

Product Specifications

CARESTREAM INSIGHT Pediatric Film
100-sheet box. Sheet Sizes – not interleaved.

Size	CAT No.
35 x 43 cm	806 4651
30 x 35 cm	177 8638
24 x 30 cm	152 4446
18 x 24 cm	190 6957
8 x 10 in.	120 6994

FOR MORE INFORMATION

Learn more about CARESTREAM INSIGHT Imaging – and what they can do for you and your patients. For other sizes or formats, contact your sales representative.