

| | | | |
|--|---|--|-------------|
| IHE Integration Statement | | Date: | 01 May 2004 |
| Vendor | Product Name | Version | |
| Eastman Kodak Company | <i>Kodak</i> RIS 2010 with Enterprise Connector | 7.0 | |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow (SWF) | DSS/Order Filler PPS Manager | None | |
| Patient Information Reconciliation (PIR) | DSS/Order Filler PPS Manager | None | |
| Presentation of Grouped Procedures (PGP) | DSS/Order Filler PPS Manager | None | |
| Internet address for Kodak IHE information: www.kodak.com/go/ihe | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| DICOM | www.kodak.com/go/dicom | | |
| Links to general information on IHE | | | |
| North America: www.rsna.org/IHE | Europe: www.ihe-europe.org | Japan: www.jira-net.or.jp/ihe-j | |