



IHE Integration Statement VNA and Vue PACS v12.1.5.5

Part # 9K2187

2016-09-02

| Integrating the Healthcare Enterprise (IHE) Integration Statement | | | |
|--|--|--|------------|
| Vendor | Product Name | Version | Date |
| Carestream Health, Inc | CARESTREAM Vue Archive (Vue for Clinical Collaboration Platform), CARESTREAM Vue PACS: Client Suite, Workflow Manager, IMS, CD Direct Suite, XDS Registry. | CARESTREAM Vue PACS, Vue Archive v12.1.5.5 | 2016-04-16 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed in this table. | | | |
| Note: This document is intended for use by Carestream customers. | | | |

| Integration Profiles Implemented | Actors Implemented | Options Implemented |
|---|--|--|
| Scheduled Workflow (SWF) | Image Manager/Image Archive Image Display PPS Manager | None |
| Imaging Object Change Management (IOCM) | Image Manager/Image Archive Change Requestor | None |
| Access to Radiology Information (ARI) | Image Manager/Image Archive Image Display Report Repository External Report Repository Access | None |
| Patient Information Reconciliation (PIR) | Image Manager/Image Archive PPS Manager | None |
| Consistent Presentation of Images (CPI) | Image Manager/Image Archive Image Display | None |
| Cross-enterprise Document Sharing (XDS.b) | Document Consumer Document Repository Document Registry Document Source | Multiple document |
| Cross-enterprise Document Sharing for Imaging (XDS.b-l) | Imaging Document Consumer Imaging Document Source | Set of DICOM Instances, Multiple documents None |

| Integration Profiles Implemented | Actors Implemented | Options Implemented |
|--|--|----------------------------|
| Cross-enterprise Document Reliable Interchange (XDR) | Document Recipient Document Source | Accept Limited Metadata |
| Cross-Enterprise Sharing of Scanned Documents (XDS-SD) | Content Consumer | None |
| Basic Patient Privacy Consents (BPPC) | Content Consumer | None |
| Patient Demographic Query HL7 V3 (PDQv3) | Patient Demographics Consumer | None |
| Patient Identifier Cross-referencing (PIX) | Patient Identifier Cross-reference Consumer | None |
| Evidence Documents (ED) | Image Display Image Manager/Image Archive | None |
| Cardiology Evidence Documents (ED-CARD) | Image Display Image Manager/Image Archive | None |
| Simple Image and Numeric Report (SINR) | Report Repository External Report Repository Access | None |
| Mammography Image (MAMMO) | Image Manager/Image Archive Image Display | None |
| Key Image Note (KIN) | Evidence Creator Image Display Image Manager/Image Archive | None |
| Nuclear Medicine Imaging (NMI) | Image Manager/Image Archive | None |
| Audit Trail and Node Authentication (ATNA) | Secured Application Secured Node | None |
| Consistent Time (CT) | Time Client | None |

| Integration Profiles Implemented | Actors Implemented | Options Implemented |
|----------------------------------|---|-------------------------------------|
| Portable Data for Imaging (PDI) | Portable Media Creator Portable Media Importer Image Display Display | Web Content None None None |

| | | |
|---|---|--|
| Link to Carestream Health IHE information: http://ihe.carestream.com | | |
| Links to Standards Conformance Statements for the Implementation | | |
| DICOM | http://dicom.carestream.com | |
| Links to general information on IHE | | |
| In North America: https://www.rsna.org/ihe.aspx | In Europe: https://www.ihe.net | In Japan: www.jira-net.or.jp/ihe-j |

Document History

| Date | Comments | Updated By |
|------------|-------------------------|---------------------|
| 2016-08-11 | Initial version | Kiran Krishnamurthy |
| 2016-09-02 | Technical writer review | Jeff Miller |

Kiran Krishnamurthy, Product Line Manager for Vue PACS – USA
For more information please contact:

| |
|--|
| ***HCIS Only*** |
| Carestream Health Technical and Innovation Center (T&IC) 150 Verona Street Rochester, NY 14608 USA |
| USA Only: 1-800-328-2910 Outside the USA: 585-627-1864 Canada: 1-866-927-1017 |
| Outside the USA and Canada: Contact your local Shared Service Center (SSC) |
| EAMER: Technical and Innovation Center (T&IC) Via al Porto Antico, Palazzina San Lorenzo 16128 - Genova -Italy In EAMER, contact the Release Management Team CNET: 215-7113 Direct Line: +39 010 2450 113 |

Authorized Representative (European Union)



Carestream Health France
1, rue Galilée
93192 NOISY-LE-GRAND CEDEX
France



Carestream



150 Verona Street
Rochester, NY 14608
USA

© Carestream Health, 2016

CARESTREAM is a trademark of Carestream Health.

Report documentation issues and feedback to HCISDocumentation@carestream.com. Include the document part number and version number (located in the footer) and any other pertinent information such as section, page, etc.