

IHE Integration Statement		Date:	October 2023
Vendor	Product Name	Version	
Carestream Health	CARESTREAM Compass X-ray System CARESTREAM DRX-1 System CARESTREAM DRX-Transportable System/Lite CARESTREAM DRX-Revolution Mobile X-ray System CARESTREAM DRX-Evolution System CARESTREAM DRX-Evolution Plus System CARESTREAM DRX-Rise System	2, Build 2.1	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
Radiation Exposure Monitoring (REM)	Acquisition Modality	None	
Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestream.com/dicom.html		
Links to general information on IHE			
North America: http://www.ihe.net	Europe: http://www.ihe-europe.net	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document
