

IHE Integration Statement			
Vendor	Product Name	Version	Date
Carestream Health	CARESTREAM PACS: CARESTREAM PACS Client, CARESTREAM Workflow Manager, CARESTREAM Information Management Solution, CARESTREAM CD Direct Suite.	V11.0	May 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display	None	
Cross-enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source Imaging Document Consumer	Set of DICOM instances, multiple documents None	
Cross-enterprise Document Sharing (XDS)	Document Repository	Multiple document	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None	
Access to Radiology Information (ARI)	Image Manager/Image Archive Report Repository Image Display Ext. Report Repository Access	None None None None	
Simple Image and Numeric Report (SINR)	Report Repository Ext. Report Repository Access	None	
Mammography Image (MAMMO)	Image Manager/Image Archive	None	
Key Image Note (KIN)	Image Manager/Image Archive	None	
Nuclear Medicine Imaging (NMI)	Image Manager/Image Archive	None	
Audit Trail and Node Authentication (ATNA)	Secured Application Image Manager/Image Archive Report Repository Ext. Report Repository Access Image Display XDS-I Document source	None	

Consistent Time (CT)	Time Client	None
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web content None None None
Link to Carestream Health IHE information: ihe.carestreamhealth.com		
Links to Standards Conformance Statements for the Implementation		
DICOM	dicom.carestreamhealth.com	
Links to general information on IHE		
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.net	In Japan: www.jira-net.or.jp/ihe-j