

| IHE Integration Statement | | Date: | October 2016 |
|--|---|---|--------------|
| Vendor | Product Name | Version | |
| Carestream Health | Carestream DRX-1 System, DRX-Mobile Retrofit Kit, DRX- Transportable, DRX-Evolution, DRX-Evolution Plus, DRX-Revolution, DRX-Ascend | 5.7 | |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow (SWF) | Acquisition Modality | None | |
| Patient Information Reconciliation (PIR) | Acquisition Modality | None | |
| Consistent Presentation of Images (CPI) | Acquisition Modality | None | |
| Portable Data for Imaging (PDI) | Portable Media Creator (PMC) | None | |
| Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| DICOM | http://www.carestream.com/dicom.html | | |
| Links to general information on IHE | | | |
| North America: http://www.rsna.org/IHE | Europe: http://www.ihe-europe.org | Japan: http://www.jira-net.or.jp/ihe-j | |

End of Document