
















## Carestream Symbols Glossary

All symbols from following standard:

ISO 15223-1 Third edition 2016-11-01 Medical devices - symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

Symbol	Symbol title	Reference Number	Meaning
	Manufacturer	5.1.1	Indicates the medical device manufacturer.
	Date of Manufacture	5.1.3	Indicates the date when the medical device was manufactured.
	Use-by date	5.1.4	Indicates the date after which the medical device is not to be used.
	Batch code	5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue number	5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Protect from heat and radioactive sources	5.3.3	Indicates a medical device that needs protection from heat and radioactive sources.
	Keep dry	5.3.4	Indicates a medical device that needs to be protected from moisture.
	Lower limit of temperature	5.3.5	Indicates the lower limit of temperature to which the medical device can be safely exposed.
	Upper limit of temperature	5.3.6	Indicates the upper limit of temperature to which the medical device can be safely exposed.
	Temperature limit	5.3.7	Indicates the temperature limits to which the medical device can be safely exposed.

	Humidity limitation	5.3.8	Indicates the range of humidity to which the medical device can be safely exposed.
	Do not re-use	5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use	5.4.3	Indicates the need for the user to consult the instructions for use.
	Caution	5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Fragile, handle with care	5.3.1	Indicates a medical device that can be broken or damaged if not handled carefully.