

Carestream

Motion Mobile X-ray System

XRS

AFFORDABLE IMAGING ON THE MOVE.

Digital-Ready. Uniquely Designed.



BRING X-RAY IMAGING RIGHT TO THE BEDSIDE – ECONOMICALLY.

When every patient who needs an imaging exam has to be transported to an X-ray room, disadvantages abound. Frail or critically ill patients can experience discomfort or trauma while being moved. Staff members can sustain injuries, especially when lifting and transporting heavy patients. Workflow slows. Costs escalate. And exams, diagnoses and treatment can all be delayed.

Now, with the Motion Mobile X-ray System, these are problems of the past. Motion Mobile is a highly affordable, analog mobile unit. It's small, lightweight and easy to use. Best of all, it allows an easy, cost-effective upgrade to digital whenever you're ready.

UNIQUELY DESIGNED WITH YOUR NEEDS – AND BUDGET – IN MIND.

Many small and mid-sized facilities have a real need for mobile X-ray capabilities. But their budgets are often stretched too tight as they are. Facilities also worry about technology obsolescence and the need to purchase new systems as their needs evolve.

The Motion Mobile System solves both of these issues. Today, its affordable purchase price makes it economically feasible for nearly any facility – and allows rapid return on investment. Tomorrow, its digital-ready design makes upgrading highly affordable, as well.

ANALOG

Choice of Manual or Anatomical Programmed Radiography (APR) mode for quick and easy selection of exam technique.

Eight-inch touch-screen monitor with easy-to-set protocols.

Pushbutton or hand-switch exposure control allows the technologists to choose their preferred method.

Tilt-step features to lift wheels over small obstacles and allow easier access to elevators.

30kW-capacitor discharge generator supports a full range of clinical applications.

Manual light-beam collimator allows fast, accurate patient positioning.

Optional DAP provides the relevant information needed for monitoring patient X-ray dose.



DIGITAL-READY.

Motion Mobile is a highly economical bridge between analog, CR and DR imaging. Its unique upgradability allows you to customize the system to meet your needs today while planning for your future. You can make a quick transition to DR whenever you're ready – and enjoy the confidence of knowing your system will have a long, useable life.

The upgraded Mobile Motion workflow is highly similar to a fully DR system – it provides fast image viewing, eliminating the longer wait time and multiple steps of using film or CR. Plus, the system automatically integrates DAP and X-ray-generator information, eliminating the need for manual entry. These efficiencies support significant cost savings.



DIGITAL

The DR upgrade kit supplements the features of our analog model with these tremendous advantages.

Hardware features include a larger, easy-to-read, 19-inch touch-screen monitor.

Carestream's DIRECTVIEW Software provides accelerated image display and easy delivery to PACS or printers.

A DAP and X-ray Generator Interface automatically acquires DAP and technique information – speeding workflow by eliminating the need for manual data entry.

Our DRX Detectors deliver a light, thin design, wireless convenience, and fast image access and cycle times.

Motion Mobile is compatible with our Family of DR Detectors

Industry-Leading Performance: Meet Our Family of DRX Detectors.

Our portfolio of DR detectors includes wireless, shareable and fixed models to meet your demands based on exam type, detector size, dosage level, and budget. Discover the right detector for your needs at carestream.com/detectors



SMART X-RAY SOLUTIONS

NO MATTER WHICH ROAD YOU'RE ON,
WE'LL GET YOU WHERE
YOU NEED TO GO.



RIGHT FOR TODAY.
READY FOR TOMORROW.

Right for Today. Ready for Tomorrow.

Carestream is ready to help you plan the most effective route to your X-ray imaging future. Our scalable equipment design and modular components mean high performance today, along with easy, affordable upgrades for years to come.

Let's plan out your equipment migration together – you'll gain the confidence that your current technology investment will continue to pay dividends well into your future.

A Community of Service and Support

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical-imaging innovation.



carestream.com

