


# DIGITAL. POWERFUL. ECONOMICAL.

Now is the Time to Step Up to DR Imaging.



# DR IMAGING MADE AFFORDABLE.

Meet the DR imaging solution designed to give private practices, clinics or any smaller imaging facility an economical way to go DR. CARESTREAM Image Suite V4 software, combined with our advanced DRX-1 System, allows you to upgrade from traditional film or CR imaging to the benefits of DR – quickly, easily and at a very attractive price.

This solution has the power to take you through the entire imaging process – from exam scheduling through archiving – with greater speed and convenience than ever before. Not to mention the spectacular image quality only Carestream DR can provide.

## CONSIDERING A CHANGE?

Upgrading to DR imaging can seem like a big step for a smaller facility. You might be wondering, “Will the change be disruptive?” or “Will it take forever for me and my staff to master the new technology?” and “Will a digital system include all of the measurement tools I need?” We can answer all your questions, and support you through a smooth and easy transition. (And, just for the record, the answers to those questions are No, No, and Yes!)

### **You'll Be Up and Running in No Time.**

Our simple installation process means you'll be fully operational before you know it. Usually in less than two hours. Plus, you can keep your existing generator and X-ray positioning equipment – no need to buy a new room.

If you're upgrading from a CARESTREAM Vita CR System with Image Suite, you'll find the DRX-1 System with Image Suite 4 features an identical user interface. So there's no need for you to learn a new one. And, the DRX-1 Detector handles just like the cassettes you're using now. This all adds up to minimal disruption and training time.

### **The CARESTREAM DRX-1 System**

The CARESTREAM DRX-1 System introduced the world's first wireless, cassette-sized detector. It protects your investment in the imaging equipment you already own – sliding right into your current system without modification, transforming it into a DR room. There's no need to modify your generator or Bucky, discard your wall stand or table, or replace your grids. You're using the equipment and interface you're already familiar with, so you'll feel right at home and be able to experience greater productivity from day one.

### **DRX Productivity**

- The wireless, battery-powered DRX Detector provides the cordless-positioning freedom of a conventional cassette.
- It quickly transmits digital images to your capture console.
- Images are rapidly ready for viewing and manipulation.
- Fast image availability can reduce exam time and recalls, improve patient satisfaction and provide quick answers.
- From the console, images can be quickly forwarded to a range of DICOM and network destinations.

### **Your Choice of DRX Detectors**

The right option for every application.

- DRX-1 Detector for general exams
- DRX-1C Detector for dose-sensitive applications
- DRX 2530C Detector with a smaller-format design for easy positioning

### **Or for Even Greater Economy, Choose a TDR Detector.**

Wired rather than wireless, Carestream's new tethered DR detectors offer the same image-processing technology, user interface and software options as our industry-leading DRX-Detectors – all at a truly affordable cost.





## IMAGE SUITE V4: THE POWER BEHIND THE PERFORMANCE.



SCHEDULE

ACQUIRE

VIEW/EDIT

DIAGNOSE/REPORT

PRINT/ARCHIVE

### Image Suite V4: a New Benchmark for Workflow Efficiency

While our "Mini-PACS" solution offers advanced image management and analysis, it's now fully optional – so if you don't need the PACS, you don't have to pay for it. Get the high-quality images you need with a system that's a superb fit, at a price that won't strain your budget.

#### Accelerated Registration

- Use Web-based software for remote-patient registration.
- Depend on our graphic menu for quick and easy selection of the body parts and views.
- Add new body parts and views, edit existing ones or apply custom processing via simple menus.

#### High-quality Images

- Fine-tune image processing with the click of a button, using simple slide bars to adjust brightness, contrast or latitude from the acquisition screen.
- Utilize specialized tools, ideal for chiropractic and orthopedic practices.
- Supports up to 2 wireless or tethered DRX-1 Systems to support the facilities X-ray needs.

#### Efficient Image Sharing and Storage

- Select optional Mini-PACS capability with a number of additional tools.
- Review images with the patient with easy tablet-viewing capability.
- Perform backups to CDs and removable USBs/networked drives for off-line storage.
- Optimize usability of images with JPEG or DICOM export tools.
- Support multifomat print composition with independent image-editing capabilities.

BOOST YOUR  
PRODUCTIVITY  
WITH THE  
THE FULL  
DIGITAL  
SOLUTION  
YOU'VE BEEN  
WAITING FOR.



**A Community of Service and Support.**

For service to depend on, look to our Customer Success Network. We work continuously to advance the performance of our imaging solutions to evolve with your changing needs – and help you make the most of your budgets and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical imaging.



[carestream.com](http://carestream.com)

