

Wireless Tablet Control.

SMART IMAGING SOLUTIONS

An Affordable Upgrade to DR – Now, with Wireless Tablet Control

Looking to upgrade your analog X-ray room or mobile unit to DR? Now you can – easily and economically – with the DRX-Transportable System/Lite.

This solution now features wireless tablet operation. There's no bulky electronics box required, and no need to work from a fixed console. Instead, a lightweight wireless tablet gives you complete control of the system and displays your images once acquired.

This allows you to move freely around the exam room, the entire radiology department or from suite to mobile units with ease.

Advantages At a Glance

- The Virtual Access Point in the detector makes this a flexible solution that enhances workflow.
- Quickly and cost effectively convert existing analog systems to DR.
- Self-powered and rechargeable, its beam-sensing technology further simplifies upgrade.
- This upgrade moves easily around the Radiology Department, from room to room or from room to mobile unit.





The Power of DR Imaging, Easier and More Portable Than Ever.

- Upgrade existing equipment to DR, instantly.
- Wireless capability lets you go anywhere.
- Capture and view images immediately at the point of patient care.
- DRX detectors are available in 3543 or 4343 wireless format with a Gadolinium (GOS) or Cesium (CSI) Scintillator.



Product Specifications

Key Components:

- Tablet PC
- i7 Processor / 12.1in touchscreen / 1.6 kg (3.6lbs)
- Rugged design to meet MIL-STD-810G and IP52 standards
- Dual Hot-Swappable battery
- DRX Detectors (sold separately)



Trusted Partner: CARESTREAM Smart Services

At Carestream, we'll be there when you need us. Our medical imaging solutions are backed by our worldwide support team of product and industry specialists – experts who work closely with you to keep your equipment up and running at peak performance. By delivering a comprehensive scope of service offerings, we provide complete peace of mind over the life of your equipment, so you can focus on what's most important – your patients.

X-Factor Thinking

The X-Factor makes the most of your capital spending with configurable, modular components powered by the portable DRX detector. Use the detector on your retrofitted mobile units on your morning rounds... and later, take it to a busy DRX-Evolution room for peak patient volumes. When the X-Factor lets you scale your solution to fit your workflow and budget, the possibilities are endless.

Industry-Leading Performance: Meet Our Family of DRX Detectors.

Our portfolio of DR detectors includes wireless, shareable and fixed models to meet your demands based on exam type, detector size, dosage, level and budget. Discover the right detector for your needs at carestream.com/detectors.

SMART IMAGING SOLUTIONS

carestream.com

