



## Consistent Quality and Affordability

Part of the Carestream Health portfolio of solutions, our Medical X-ray Cassettes with Green 400 Screens are designed for medical x-ray green-sensitive films. These films are orthochromatic and suitable for all types of radiographic examinations. And they provide consistent quality in an affordable cassette for general radiography.

## ADVANTAGES AT A GLANCE

- Consistent performance
- Multiple sizes
- Sharp images
- Durable design



# Get Predictable Results From a Name You Trust.

## The Right Cassette for The Job

Our high-quality cassettes are designed to suit the way you work. With 9 sizes to choose from, you can be sure to find the right-sized cassette for every examination.

## 400-Speed Rare-Earth Green Screens

Our proven medical x-ray cassettes feature fast, 400-speed green-emitting screens for sharp images and low patient exposure to radiation.

## Tight, Uniform Film Contact

A high-quality back panel and specially fitted foam provide tight, uniform contact between the film and screen. This ensures sharp, high-quality images for every examination.

## With or Without Window ID

Choose the x-ray cassette that best meets your needs—available with or without a window for patient identification.

## Built-in Backscatter Control

The lead construction on the back panel helps control backscatter radiation and contributes to sharp, clear diagnostic images.

## Durability

These affordable cassettes are designed for demanding imaging environments. Aluminum alloys are used for the front plate and back panel, providing rigidity. Reinforced corners add strength. Cassettes are also scratch-resistant and covered by a limited one-year warranty.

## Easy Loading and Unloading

A simple, recessed latch design makes it easy for you to load and unload film after film.



## CATALOG NUMBERS

Size	With Window	Without Window
8 x 10 inch	6805642	6807093
10 x 12 inch	6805659	6807101
11 x 14 inch	6805667	6807127
18 x 24 cm	6805709	6807002
18 x 43 cm	6805717	6807010
24 x 24 cm	6805725	6807028
24 x 30 cm	6805733	6807036
35 x 35 cm	6805766	6807069
35 x 43 cm	6805774	6807077



[carestream.com](http://carestream.com)

