



Reliability And Value For Mammography.

The CARESTREAM Mammography X-Ray Processor offers reliable automated processing at a very affordable price. This tabletop unit fits almost anywhere in your clinic, private practice, or other health facility. A simple, streamlined design provides high reliability, ease of use, and easy service maintenance.

This Mammography X-Ray Processor delivers consistent performance at attractive prices. It's all part of Carestream Health's commitment to help every customer meet their image quality, productivity, and service needs.

CARESTREAM Mammography X-Ray Processor (MXP-M)

Image Quality

Designed for mammography applications, this tabletop unit can deliver reliable image quality for a shallow-tank processor – especially when used with Carestream medical x-ray processing chemicals.

Consistent Results

Reduce artifacts and obtain consistent results with the aid of anti-crystallization and anti-oxidation cycles that help maintain clean racks and proper chemical activity when the processor is not in use.

Films Up To 35cm Wide

Process sheet films up to 35 cm x 43 cm or as small as 18 cm x 24 cm.

Quick Turnaround

View images in as little as 105 seconds. Compare this to 5 minutes for standard manual processing.

Chemical, Water, And Energy Savings

Reduce chemical use and cost with semi-automatic replenishment that uses a microsensor to replenish chemicals on the basis of actual film length processed. Save water and energy costs with an integrated standby mode that reduces consumption when the processor is not in use.

Easy To Use And Maintain

A convenient control panel enables easy filling of tanks during the warm-up cycle. Temperature setting is simple and accurate with microprocessor control.

The processor automatically signals when it is ready to process the next film. Few moving parts and easy-to-service film-transport racks simplify maintenance and service.

Service And Support

Authorized distributors of Carestream products who sell the CARESTREAM Mammography X-Ray Processor provide service and support.

Streamlined Installation

Ease of installation means you'll be processing films quickly. Facility staff can install the processor, or you can request that a trained Carestream Health-authorized distributor complete the installation.

Choice Of Optional Accessories

Enhance the convenience and value of your Mammography X-Ray Processor with a choice of optional accessories:

- Medical X-Ray Processor Stand
CAT No. 870 0338
- Medical X-Ray Processor Cabinet
CAT No. 184 3226
- Two Replenishment Tanks, 25 L
CAT No. 855 4313

Product Specifications

- **Processing Cycle**
105 second cycle capability
- **Throughput**
89 films/hour (18 cm x 24 cm)
71 films/hour (24 cm x 30 cm)
@ 105-second cycle
- **Dimensions**
Length: 77 cm/30 in.
Width: 59 cm/23 in.
Height: 42 cm/16.5 in.
- **Weight**
Without solutions: 35 kg
With solutions: 50 kg
- **Tank volumes**
Developer: 5.0 L
Fixer: 5.0 L
Wash: 5.0 L
Replenishment tanks, 12 L included
- **Electrical**
Model 112: 220V - 240V, 50Hz, 16 Amp
Model 113: 220V - 240V, 60Hz, 15 Amp
Conforms to CE standards
Consumption:
- Processing: 1.4 kW/hr (average),
1.99 kW/hr (maximum)
- Standby: 0.12 kW/hr
- **Catalog numbers**
Model 112:
CAT No. 516 1336 - 230V 50Hz
Model 113:
CAT No. 516 1344 - 230V 60 Hz

www.carestream.com



"Rx only"

© Carestream Health, Inc., 2018. CARESTREAM is a trademark of Carestream Health. CAT 864 5145 12/18
Mammography X-Ray Processor (MXP-M) not available in all countries.



Carestream