



Exceptional Quality and Value

Controlling costs is imperative. So is quality care. Carestream understands you want both — affordability and the tools you need to provide consistent, quality care. This understanding has led to the introduction of Select Cassettes. Part of the Carestream Health portfolio of solutions, Select Cassettes with LANEX Screens are designed for seamless compatibility with green-sensitive medical x-ray films.

ADVANTAGES AT A GLANCE

- Exceptional value
- Same great screens
- Clean, sharp images
- One-year warranty



CARESTREAM Select Cassettes with LANEX Screens

Make No Compromise.

High-Quality Cassettes at an Affordable Price

If you are looking for the perfect mix of value and performance, new Select Cassettes deliver. These cassettes feature high-quality construction, based on years of success in healthcare facilities around the world. They deliver the same great performance characteristics you expect from a leader in imaging technology.

Superb Image Quality

Select Cassettes feature our best-selling LANEX Screens in ten common sizes. As always, you can count on our proven green-emitting screens for highly consistent and exceptionally sharp images. You can choose from our line of regular, fine or fast screens. In addition to superb image quality, our innovative screen design helps reduce patient exposure to radiation.

Uniform Screen-Film Contact

Select Cassettes feature a specially designed back panel and fitted foam lining to provide uniform contact between the film and screen. You can count on high-quality images exam after exam.

Quality Construction

The back panel of the new Select Cassettes features the lead lining construction that minimizes backscatter. This proven design contributes to sharp, clear diagnostic images.



CATALOG NUMBERS

Size	Regular Screens (with window)	Regular Screens (without window)	Fine Screens (with window)	Fast Screens (with window)
8 x 10 inch	1005172	1017847		1005438
10 x 12 inch	1004761	1017870		1005222
11 x 14 inch	1004829	1017896		
18 x 24 cm	1004828	1017672		1005263
18 x 43 cm	1004977	1017813		1005289
24 x 30 cm	1005016	1017706	1005479	1005339
30 x 35 cm	1005032	1017722		1005354
30 x 40 cm	1005073	1017797		1005362
35 x 35 cm	1005081	1017748		1005404
35 x 43 cm	1005156	1017763		1005412



carestream.com

