



Big Performance. Small Package.

The CARESTREAM DRX-1 System unleashed the power of the X-Factor, with the world's first wireless, cassette-sized detector. Its advanced platform allowed a single DRX detector to work seamlessly across the entire family of DRX products.

Soon after, we introduced the second-generation DRX detector, a lighter, faster detector with enhanced features.

Now the next star in the DRX family is close at hand. The DRX 2530C Detector* will retain all the high-performance benefits of our current detector in a smaller-format design – optimized for pediatric and orthopedic imaging.

*Not Commercially Available

ADVANTAGES AT A GLANCE

The DRX 2530C Detector offers:

- Fast, easy positioning in pediatric incubator trays
- Cesium iodide (CsI) design for dose-sensitive pediatric applications
- Ideal for orthopedic tabletop imaging, due to its flexible positioning and easy handling
- Sharing across the entire DRX family of products
- Proven technology, with over 5,000 detectors in use world wide



CARESTREAM DRX 2530C Detector



DRX 2530C Specifications:

Wireless standard 802.11n

Image size: 25 x 30 cm

Pixel pitch: 0.139 mm

External dimensions:

28.4 x 41.5 x 1.6 cm (11.2 x 16.3 x .63 in)

Additional Detector Benefits:

- Wireless performance, to eliminate the hassles and hazards of cables and reduce the risk of infection
- Fast performance and instant access to images
- Ergonomically designed handle for fast and safe positioning
- Compatible with the DRX-1 Detector battery charger to help eliminate additional equipment
- Fits into the DRX-Revolution small-format detector holder

Protection Against the Unexpected

The DRX 2530C will be eligible for coverage under the DRX Accident Protection Plan. This uniquely designed plan can be added to your standard Carestream DRX Warranty or Service Agreement. The plan will allow replacement of your DRX 2530 Detector at a fraction of the full cost, should damage occur due to an accident in the course of normal clinical use. For you, this means added confidence and peace of mind, convenient, one-call resolution and fast replacement to maximize uptime.

Meet The Family

The imaging solutions in the DRX Family are all designed around the remarkable X-Factor wireless detector platform. This allows the detector to be shifted from any DRX system to any other one – and also work with your existing equipment.

- **The DRX-1 System** easily converts existing analog rooms to DR.
- **The DRX-Evolution** is a full DR suite configurable to meet your workflow, budget and space.
- **The DRX-Ascend**, a floor-mounted DR room, provides advanced capabilities at a very affordable price.
- **The DRX-Revolution** is our highly maneuverable, compact mobile system packed with innovative features, for wireless DR imaging wherever it's needed.
- **The DRX-Mobile Retrofit Kit** lets you upgrade to DR capability while still leveraging your investment in your current mobile systems from GE, Siemens and Shimadzu.
- **The DRX-Transportable/Universal Mobile** makes it easy to upgrade mobile systems from many other manufacturers.
- **The DRX-Transportable/Field Portable** allows you to take reliable DR performance into the field. It's an all-in-one solution, packed in a rugged case.



carestream.com/drx

