

IHE Integration Statement		Date:July 2013
Vendor	Product Name	Version
Carestream Health Inc.	<i>DryView 7920 Laser Imager</i>	1.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles implemented	Actors Implemented	Options Implemented
Consistent Presentation of Images (CPI)	Print Server	None
Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html		
Links to Standards Conformance Statements for the Implementation		
DICOM	http://www.carestream.com/dicom.html	
Links to general information on IHE		
North America: http://www.rsna.org/ihe.aspx	Europe: http://www.ihe-europe.net	Japan: http://www.jira-net.or.jp/ihe-j

End of Document