

Bring the key capabilities of a radiology room right to the patient – at the bedside, in the OR, the ICU or EDU. For maximum ease and convenience, be sure to put all of these capabilities to work.



## CARD/BADGE READER

Quickly power up and start your rounds with a single badge swipe. Many pcProx(R) 125 kHz cards and AIR ID(R) 13.56 Mhz cards are supported free of charge as a standard feature. Optional proximity badges may be ordered by customers not using a standard supported badge.

## PATIENT DATA PRIVACY FEATURE

With the press of a button, you can keep confidential patient information from showing on the column screen.

## LOCKING MECHANISM FOR DETECTORS

Protect your investment and use the locking feature when detectors are not in use. This will ensure they stay with the mobile unit and are not misplaced.

## DETECTOR HOLDER FOR BAGGING DETECTORS

Utilize the detector holder to quickly and easily bag detectors prior to imaging. This protects both the patient and the detector from bacterial contamination.

## LOCKABLE STORAGE BINS

Pack extra batteries, paperwork, markers, bags and gloves to have on hand, whenever needed.



## GRID ALIGNMENT SYSTEM

This unique feature makes tube-and-grid alignment easier than ever before – helping to deliver superb X-ray quality.

## BUMPER BAR

Remember, this feature is here to protect the unit from damage in case of collision.

## UNIT-CHARGING POWER CORD

Plug in at night so you're always charged and ready to go in the morning. To prevent damage to the power cord and extend its life, be sure to position the side of the unit with the power cord and plug close to the wall outlet.



## CARESTREAM SERVICE

Count on our team of service professionals and service offerings to give you peace of mind and let you focus on what's most important – your patients.