



Read Digital Breast Tomosynthesis from your Vue Mammo Workstation.

Whether you're reading traditional Mammograms, Breast Ultrasound, Breast MRI, General Radiology or Digital Breast Tomosynthesis (DBT) exams, the CARESTREAM Vue Mammo Workstation lets you do it all from a single desktop. Carestream's DBT module adds unique workflow capabilities and specialized tools that optimize the reading of Digital Breast Tomosynthesis exams at the same time as other procedures.

In addition, Carestream can securely archive your Digital Breast Tomosynthesis exams for accessibility from any Vue Mammo Workstation. Choices include enterprise storage with CARESTREAM Vue Archive or cloud storage with CARESTREAM Vue Cloud.

ADVANTAGES AT A GLANCE

- Read all procedures from a single desktop.
- The workstation is optimized for Digital Breast Tomosynthesis.
- Easily add this optional module to a new or existing Vue Mammo Workstation or Vue PACS.
- Secure archive solutions make your procedures available from any Vue Mammo Workstation.

CARESTREAM Vue Mammo Workstation Digital Breast Tomosynthesis Module

Streamlined Workflow

Store, route and display Digital Breast Tomosynthesis exams using industry-recognized standards.

- Storage of Digital Breast Tomosynthesis exams from DICOM-compliant acquisition devices
- Query / Retrieve Digital Breast Tomosynthesis exams by DICOM-compliant devices
- Worklist Indicator for quick identification of exams with Digital Breast Tomosynthesis images

Optimized Comparison

Compare multiple prior exams and move images as you read.

- User-configurable Hanging Protocols to display Digital Breast Tomosynthesis procedures along with any other procedure
- Entire exam can be replaced with a single selection
- Drag and Drop images as you read with automatic resizing and positioning

Specialized Tools

Save time while enhancing diagnostic confidence.


- Automatic positioning of Digital Breast Tomosynthesis images and mammograms to help eliminate manual manipulations
- Automatic Same Sizing of Digital Breast Tomosynthesis and mammograms for ease of comparing a change in pathology
- Automatic skin-line detection to focus operations to the tissue area
- Digital Breast Tomosynthesis exams viewable in cine mode or through self-paced image scrolling
- Concurrent Magnifying Glasses provide a close-up comparison of pathology across multiple views and procedures
- Cross-reference lines to aid in locating suspicious areas on opposing views
- Save and catalog annotations and key images

Trusted Partner

Carestream is your trusted partner, providing a suite of products and solutions to support your mammography needs. Look to Carestream to help you provide superb patient care.

More Information

To learn more about Carestream Digital Mammography Solutions and our partners, contact your local Carestream representative or authorized dealer.



Digital Breast Tomosynthesis

2D Mammogram



carestream.com/mammo

