



Streamline Workflow With Linked CR and DR Systems

The DIRECTVIEW Capture Link System links up to five DIRECTVIEW DR and/or CR Systems to improve your workflow and increase productivity in high-volume areas. Now you can integrate DR and CR systems together for greater flexibility when a study includes both DR and CR images. Capture Link makes sharing cassette identification, scanning and performing image review functions on linked CR systems easier for all.

Add up to ten DIRECTVIEW Remote Operations Panels (ROPs) or PCs with DIRECTVIEW Remote Patient Data Entry Software (RPDES) for even greater convenience.



DRX-Evolution System


DIRECTVIEW Elite CR System

ADVANTAGES AT A GLANCE

- Identify, scan, and review CR images anywhere
- Scan cassettes quickly at any linked DIRECTVIEW CR system
- Cassettes from one study can be scanned at one DIRECTVIEW CR system or at multiple systems*
- Review, check quality and deliver images to final destinations from any linked DIRECTVIEW CR system monitor or ROP.
- Identify cassettes at up to five linked DIRECTVIEW CR systems—and up to ten wall-mounted DIRECTVIEW ROPs

* Cassettes from a long-length image must be scanned on the same CR system

CARESTREAM DIRECTVIEW Capture Link System



Managing studies with DR and CR images has never been easier

- Link DIRECTVIEW DR and CR systems together seamlessly
- Identify CR cassettes at the DR console and scan them at a remote or in-room CR reader
- View the entire study at the DR console—DR and CR images—at once for greater productivity and convenience
- Reprocess, approve and redirect CR images from your DIRECTVIEW DR system operator console
- Forward DR and CR images together to the final output destinations

Harness the power of multiple systems

- Flexible configuration fits the way you work
- Support better workflow and increased throughput in mid- to high-volume, centralized radiography applications
- Modify linked systems to meet changing workflow needs—add or remove DR and CR systems, ROPs or RPDES connections
- DIRECTVIEW DR and CR systems share a consistent, intuitive user interface to minimize training

Improve CR data sharing and access

- Access worklists and all patient and exam information from any linked DIRECTVIEW CR system or device
- Update a single patient record, reprocess images, check for quality and deliver to the output destinations from any linked CR system or device

DIRECTVIEW Capture Link System requirements

- Access worklists and all patient and exam information from any linked DIRECTVIEW CR system or device
- Update a single patient record, reprocess images, check for quality and deliver to the output destinations from any linked CR system or device

Simple Capture Link:

- Links two CRs or a CR and DR systems
- With both systems running DIRECTVIEW software, one system is selected to be the server for the two units
- Includes all the functions of DIRECTVIEW Capture Link
- Allows addition of up to ten Remote Operations Panels

