

<b>IHE Integration Statement</b>		<b>Date:</b>	15 Oct 2004
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Eastman Kodak Company	<i>Kodak DirectView PACS System 5</i>	5.2	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display Print Composer	None	
Access to Radiology Information (ARI)	Image Manager/Image Archive Report Repository Image Display	None	
Simple Image and Numeric Report (SINR)	Report Repository	None	
Key Image Note (KIN)	Image Manager/Image Archive	None	
Basic Security (SEC)	Secure Node Image Manager/Image Archive Image Display	None	
<b>Internet address for Kodak IHE information:</b> <a href="http://www.kodak.com/go/ihe">www.kodak.com/go/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>DICOM</b>	<a href="http://www.kodak.com/go/dicom">www.kodak.com/go/dicom</a>		
<b>Links to general information on IHE</b>			
North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	