

IHE Integration Statement		Date:	February 2018
Vendor	Product Name	Version	
Carestream Health Inc.	MyVue Center	1.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Composer Print Server	None	
Access to Radiology Information (ARI)	Image Display	None	
Portable Data for Imaging (PDI)	Portable Media Creator	USB Media	
Consistent Time (CT)	Time Client	None	
Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestream.com/dicom.html		
Links to general information on IHE			
North America: http://www.rsna.org/ihe.aspx	Europe: http://www.ihe-europe.net	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document

AC7050