



## AFFORDABLE IMAGING ON THE MOVE.

When every patient who needs an imaging exam has to be transported to an X-ray room, disadvantages abound. Frail or critically ill patients can experience discomfort or trauma while being moved. Workflow slows. Costs escalate. And exams, diagnoses and treatment can all be delayed.

Now, with the Motion Mobile X-ray System, these are problems of the past. Motion Mobile is a highly affordable, analog mobile unit. It will change the way you think about mobile imaging.



## ADVANTAGES AT A GLANCE

- Light, compact and easy to transport.
- Intuitive touch-screen for ease of operation.
- Rapid and precise positioning, thanks to a manual light-beam collimator.
- Convenient storage bin keeps cassettes close at hand.

# Motion Mobile X-ray System

## UNIQUELY DESIGNED WITH YOUR NEEDS – AND BUDGET – IN MIND.

Many small and mid-sized facilities have a real need for mobile X-ray capabilities. But their budgets are often stretched too tight as they are. The Motion Mobile System solves both of these issues, making it economically feasible for nearly any facility. Consider this advanced feature set – and consider the value.



Manual light-beam collimator allows fast, accurate patient positioning.

Optional DAP provides the relevant information needed for monitoring patient X-ray dose.



Choice of Manual or Anatomical Programmed Radiography (APR) mode for quick and easy selection of exam technique.

Eight-inch touch-screen monitor with easy-to-set protocols is simple to learn and simple to use.

30kW-capacitor discharge generator supports a full range of clinical applications.

Tilt-step features to lift wheels over small obstacles and allow easier access to elevators.

[carestream.com](http://carestream.com)

