

IHE Integration Statement		Date:	Dec 2011
Vendor	Product Name	Version	
Carestream Health	Carestream <i>DirectView</i> DR Systems (DR 9500, DR 3500, DR 7500) Carestream DRX-1 System, DRX-Evolution, DRX-Mobile Retrofit Kit,	5.5	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
Portable Data for Imaging (PDI)	Portable Media Creator (PMC)	None	
Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestream.com/dicom.html		
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe: http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document