

ClinAsept



Carestream Dental Films with a ClinAsept Barrier

INSIGHT and Ultra-speed films with a ClinAsept Barrier have a built-in barrier to help manage infection control. Each intraoral packet is encased in a clear barrier that guards the film packet against contaminants that may be present in oral fluids. Use these products to minimize the risk of cross-contamination in your office.




Minimize Cross-Contamination and Maximize Confidence

For real time savings and peace of mind, use INSIGHT or Ultra-speed film packets with ClinAsept barriers. Available in single- or double-film format, for size 2 and single-film format for size 0. Each packet comes safely prepackaged with its own barrier. You don't have to do anything extra before taking the radiograph, except feel confident that the risk of cross-contamination in the darkroom is minimized.

Get in Touch Today

Make INSIGHT or Ultra-speed film packets with ClinAsept barriers part of your key infection control practices. For speed and convenience, they're a real plus in your dental practice.

How to use prepackaged films with barriers



Ordering Information

Code	Size	Quantity Per Package	REF No
INSIGHT dental film with a ClinAsept Barrier			
IP-21C	Size 2	100 1-film packets	153 9931
IP-22C	Size 2	100 2-film packets	156 0390
IP-01C	Size 0	75 1-film packets	120 0328
Ultra-speed dental film with a ClinAsept Barrier			
DF-58C	Size 2	100 1-film packets	171 7131
DF-57C	Size 2	100 2-film packets	898 7877
DF-54C	Size 0	75 1-film packets	144 5360

For more information, call **800.933.8031** or visit **carestream.com**.