

<b>IHE Integration Statement</b>		<b>Date:</b>	Apr 2021
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Carestream Health	CARESTREAM DRX-Ascend System Quantum Q-Rad-DIGITAL DRX-Series CARESTREAM Compass X-ray System CARESTREAM DRX-1 System CARESTREAM Mobile X-ray System CARESTREAM DRX-Transportable System/Lite CARESTREAM DRX-Revolution Mobile X-ray System CARESTREAM DRX-Revolution Nano Mobile X-ray System CARESTREAM DRX-Evolution System CARESTREAM DRX-Evolution Plus System	1, Build 1.9	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
Radiation Exposure Monitoring (REM)	Acquisition Modality	None	
<b>Internet address for Carestream Health IHE information:</b> <a href="http://www.carestream.com/ihe.html">http://www.carestream.com/ihe.html</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>DICOM</b>	<a href="http://www.carestream.com/dicom.html">http://www.carestream.com/dicom.html</a>		
<b>Links to general information on IHE</b>			
North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	

End of Document

---