


Carestream

Specialized
Equine Imaging



DIGITAL. POWERFUL. CUSTOMIZED.

Solutions Designed with
a Vet's Needs in Mind.



A FULLY DIGITAL SOLUTION, SPECIALIZED FOR VETS.

In the barn or on the go, equine imaging has its own special requirements. Here's everything you need to leverage the benefits of portable, high-quality, digital X-ray.

DRX-Transportable System / Field Portable

A full digital X-ray system packed in a rugged, water- and dust-resistant case:

- Ideal for in the field, in the barn, and for easy transport
- Compatible with mobile X-ray equipment from multiple suppliers
- Self-powered and rechargeable
- Uses CARESTREAM DRX Detectors

DRX Detectors

Carestream offers a range of detector options, so you can select the ideal balance of performance and price for your veterinary clinic. All are available with your choice of either a Gadolinium (GOS) or a Cesium (CsI) scintillator for dose-sensitive applications.





CARESTREAM IMAGE SUITE 4 SOFTWARE.

Many digital-imaging systems marketed to vets come with a generic, human-based clinical interface that you need to adapt. By contrast, Image Suite software offers an interface designed for the way vets work. And the entire system is loaded with custom features you're looking for.

- Our order-entry screens have custom fields for selection of species.
- You can input unique identifiers like tattoos and chip IDs.
- Reference images enable accurate positioning, while sample X-rays provide guides for enhanced image quality.
- Our animal-image database fine-tunes processing parameters for specific species, body parts and views.
- Our Laminitis Measurement Tool helps you arrive at a precise diagnosis faster than ever by measuring the rotation angle of the pedal bone relative to the dorsal hoof wall. This view gives useful information about the degree of pedal rotation, sole thickness, dorsal hoof wall thickness and distal displacement of the pedal bone.

Industry-Leading Performance: Meet Our Family of DRX Detectors.

Our portfolio of DR detectors includes wireless, shareable and fixed models to meet your demands based on exam type, detector size, dosage level, and budget. Discover the right detector for your needs at carestream.com/detectors

Faster Results

This solution delivers the high-quality images you need in minutes or less – so diagnoses can be made far more quickly. Plus, reviewing diagnostic images on-screen with animal owners makes it far easier to explain the animal's condition with clarity and speed – and show owners why your recommended treatment is the best option. Of course, rapid image access also means expedited workflow and enhanced productivity in your practice.

Ease of Use

The DRX-Transportable System with Image Suite 4 is exceptionally simple to learn – training time is minimal, and the workflow matches the way you handle film today. Moreover, the system's interface is highly intuitive, so you and your staff will be using it like experts in no time.

Imaging You Can Depend On

You'll benefit from sharp, clear images that support more accurate review and diagnosis. They've been fine-tuned for the body parts and views of the animals you treat. So don't settle for processing designed for a human ankle – when what you really need is processing optimized for a stifle.



BOOST YOUR
PRODUCTIVITY
WITH THE
DIGITAL
SOLUTION
VETS HAVE BEEN
WAITING FOR.



A Community of Service and Support

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical-imaging innovation.



carestream.com

