

<b>IHE Integration Statement</b>		<b>Date:</b>	Jun 2005
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Eastman Kodak Company	<i>Kodak</i> RIS 2010 Software with Enterprise Connector	8.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	DSS/Order Filler PPS Manager	Appointment Notification	
Patient Information Reconciliation (PIR)	DSS/Order Filler PPS Manager	None	
Basic Security (SEC)	Secure Node DSS/Order Filler PPS Manager	None	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-referencing Consumer	None	
Retrieve Information for Display (RID)	Information Source Display	None	
<b>Internet address for Kodak IHE information:</b> <a href="http://www.kodak.com/go/ihe">www.kodak.com/go/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.hl7.org">www.hl7.org</a>		
<b>DICOM</b>	<a href="http://www.kodak.com/go/dicom">www.kodak.com/go/dicom</a>		
<b>Links to general information on IHE</b>			
North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	