

Carestream

Motion Mobile X-Ray System

# Mobile. Affordable. Upgradable.

Uniquely Designed. Digital Ready.



SMART IMAGING SOLUTIONS

# Affordable Imaging on the Move.

When every patient who needs an imaging exam must be transported to an X-ray room, disadvantages abound. Frail or critically ill patients can experience discomfort or trauma while being moved. Meanwhile, workflow slows. Costs escalate. And exams, diagnoses and treatment can all be delayed.

Now, meet the Motion Mobile X-ray System. It's a highly affordable, analog mobile unit – small, lightweight and easy to use.

Best of all, it allows an easy, cost-effective upgrade to digital whenever you're ready – with the simple addition of our DRX-Transportable/Lite. It uses a lightweight, wireless tablet to give you complete control of the system and displays your images once you acquire them.

## Advantages at a Glance

- Light, compact and easy to transport.
- Intuitive touchscreen for ease of operation.
- Rapid, precise positioning, thanks to a manual light-beam collimator.
- Fast, easy upgrade to digital with the addition of our lightweight, wireless tablet.



# Designed with Your Needs – and Budget – in Mind.

Many small- and mid-sized facilities have a real need for mobile X-ray capabilities. But their budgets are often stretched too tight as they are. The Motion Mobile System solves both of these issues, making it economically feasible for nearly any facility. Consider its advanced feature set – and consider the value.

## Go Lite!

### Go Digital, Whenever You're Ready.

Simply add our DRX-Transportable System/Lite and its wireless tablet, with a DRX Plus or DRX Core Detector, and Motion Mobile becomes a full digital system – with the image quality and accelerated workflow only DR imaging can provide. This makes Motion Mobile a highly economical bridge from analog or CR to the performance power of DR.

Bottom line? You can make a quick transition to DR at any time – with no fear of technology obsolescence! You'll enjoy the confidence of knowing your system will have a long, useable life well into the future.



The DRX-Transportable System/Lite features wireless tablet operation. There's no bulky electronics box required, and no need to work from a fixed console – you can move freely around the exam room. A DRX-Detector and battery charger complete the system.

# Feature Rich. Digital Ready.



Manual light-beam collimator allows fast, accurate patient positioning.

Optional DAP provides the relevant information needed for monitoring patient X-ray dose.

Choice of Manual or Anatomical Programmed Radiography (APR) mode for quick and easy selection of exam technique.

Eight-inch touchscreen monitor with easy-to-set protocols is simple to learn and use.

Tilt-step feature to lift wheels over small obstacles and allow easier access to elevators.

30kW-capacitor discharge generator supports a full range of clinical applications.



## Tap the Power of X-Factor Thinking.

Our wireless, lightweight DRX Detector not only controls Motion Mobile – it's also shareable with our entire line of DRX X-ray rooms and mobile systems. And, it features a Virtual Access Point that enhances workflow.

NO MATTER WHICH ROAD YOU'RE ON,  
WE'LL GET YOU WHERE  
YOU NEED TO GO.

SMART IMAGING SOLUTIONS

### Right for Today. Ready for Tomorrow.

Carestream is ready to help you plan the most effective route to your X-ray imaging future. Our scalable equipment design and modular components mean high performance today, along with easy, affordable upgrades for years to come. Let's plan out your equipment migration together – you'll gain the confidence that your current technology investment will continue to pay dividends well into your future.

### A Community of Service and Support

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical-imaging innovation.



[carestream.com/dr](http://carestream.com/dr)

"Rx only"



© Carestream Health, Inc., 2017. CARESTREAM and DIRECTVIEW are trademarks of Carestream Health. CAT 9K4019 09/17

**Carestream**