

Wireless Tablet Control.

SMART IMAGING SOLUTIONS

Upgrade to Affordable DR with an All-In-One Solution.

Looking to upgrade your analog X-ray room or mobile unit to DR? Now you can – easily and economically – with the CARESTREAM DRX-Transportable System/Lite.

This solution operates with a versatile, wireless tablet. There's no bulky electronics box required, and no need to work from a fixed console. Instead, a lightweight, wireless tablet gives you complete control of the system and displays your images once acquired.

The DRX-Transportable System/Lite allows you to move freely around the exam room, the entire radiology department or from suite to mobile units with ease.

Advantages At a Glance

- Provides a lightweight, portable solution with no need for external components.
- Quickly and cost-effectively converts existing room and mobile analog systems to DR.
- Simplifies your upgrade because it's self-powered and rechargeable, with beam-sensing technology.
- Moves easily around the Radiology Department, from room to room or from room to mobile unit.





The Power of DR Imaging, Easier and More Portable Than Ever.

- Upgrade existing equipment to DR, instantly.
- Wireless capability lets you go anywhere.
- Capture and view images immediately at the point of patient care.
- Choose from a variety of CARESTREAM DR Detector options:

DRX Plus These versatile detectors feature the X-Factor and are available in multiple sizes, including a 25x30 small-format detector that's ideal for pediatric imaging. **Lux 35** Our uniquely designed glass-free detector includes ergonomic features for easier positioning and more comfortable patient exams.

Focus 35C and 43C These provide economical ways to upgrade to full wireless DR technology.



Product Specifications

Key Components:

- Tablet PC: HP Elite x2 G4 Tablet
 - 8th Generation Intel® Core™ i5 Processor
 - 13.0 inch Multi-Touch Display with 3000 x 2000 resolution
 - 16GB of Memory
 - 512GB Solid-State Drive for Storage
 - 1.4 kg (3.1 lbs)
- DR Detectors (sold separately)

Trusted Partner: CARESTREAM Smart Services.

At Carestream, we'll be there when you need us. Our medical imaging solutions are backed by our worldwide support team of product and industry specialists – experts who work closely with you to keep your equipment up and running at peak performance. By delivering a comprehensive scope of service offerings, we provide complete peace of mind over the life of your equipment, so you can focus on what's most important – your patients.

X-Factor Thinking.

The X-Factor makes the most of your capital spending with configurable, modular components powered by the portable DRX and Lux detectors. Use the detector on your retrofitted mobile units on your morning rounds...and later, take it to a busy DRX-Evolution room for peak patient volumes. When the X-Factor lets you scale your solution to fit your workflow and budget, the possibilities are endless.

Industry-Leading Performance:

Meet Our Family of DR Detectors.

Our portfolio of DR detectors includes wireless, shareable, glass-free and fixed models to meet your demands based on exam type, detector size, dosage level and budget. Discover the right detector for your needs at carestream.com/detectors.



carestream.com



"Rx only"

© Carestream Health, Inc., 2022. CARESTREAM is a trademark of Carestream Health. CAT 2000 214 02/22

