

IHE Integration Statement		Date:	Dec 2019
Vendor	Product Name	Version	
Carestream Health	CARESTREAM DRX-Ascend System Quantum Q-Rad-DIGITAL DRX-Series	1.0, Build 1.5	
Carestream Health	CARESTREAM Compass OTC system	1, Build 1.6	
Carestream Health	CARESTREAM DRX-Ascend System Quantum Q-Rad-DIGITAL DRX-Series CARESTREAM DRX-1 System CARESTREAM Mobile X-ray System CARESTREAM DRX-Revolution Mobile X-ray System CARESTREAM DRX-Revolution Nano Mobile X-ray System CARESTREAM Motion Mobile X-ray System CARESTREAM DRX-Evolution System CARESTREAM DRX-Evolution Plus System	1, Build 1.7	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
Radiation Exposure Monitoring (REM)	Acquisition Modality	None	
Internet address for Carestream Health IHE information: http://www.carestream.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestream.com/dicom.html		
Links to general information on IHE			
North America: http://www.ihe.net	Europe: http://www.ihe-europe.net	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document
