

| <b>IHE Integration Statement</b>   |  |   |             |
|--|--|---|-------------|
| <b>Vendor</b>  | <b>Product Name</b>  | <b>Version</b>  | <b>Date</b> |
| Carestream health  | Carestream PACS: <a href="#">Client Suite, Workflow Manager and CD Direct Suite, IMS.</a>  | CSP v10.2   | March 2008  |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |  |   |             |
| <b>Integration Profiles Implemented</b>  | <b>Actors Implemented</b>  | <b>Options Implemented</b>                            |             |
| Scheduled Workflow (SWF)   | Image Manager/Image Archive<br>Image Display<br>PPS Manager  | None  |             |
| Patient Information Reconciliation (PIR)   | Image Manager/Image Archive<br>PPS Manager   | None  |             |
| Consistent Presentation of Images (CPI)  | Image Manager/Image Archive<br>Image Display   | None  |             |
| Cross-enterprise Document Sharing for Imaging (XDS-I)  | Imaging Document Source<br><br>Imaging Document Consumer   | Set of DICOM Instances,<br>Multiple documents<br>None |             |
| Cross-enterprise Document Sharing (XDS)  | Document Repository  | Multiple document                                     |             |
| Patient Identifier Cross-referencing (PIX)   | Patient Identifier Cross-reference Consumer  | None  |             |
| Access to Radiology Information (ARI)  | Image Manager/Image Archive<br>Report Repository<br>Image Display<br>Ext. Report Repository Access   | None<br>None<br>None<br>None                          |             |
| Simple Image and Numeric Report (SINR)   | Report Repository<br>Ext. Report Repository Access   | None  |             |
| Mammography Image (MAMMO)  | Image Manager/Image Archive  | None<br>None  |             |
| Key Image Note (KIN)   | Image Manager/Image Archive  | None  |             |
| Nuclear Medicine Imaging (NMI)   | Image Manager/Image Archive  | None  |             |
| Audit Trail and Node Authentication (ATNA)   | Secured Application<br>Image Manager/Image Archive<br>Report Repository<br>Ext. Report Repository Access<br>Image Display<br>XDS-I Document source | None  |             |
| Consistent Time (CT)   | Time Client  | None  |             |

|  |   |  |
|--|---|--|
| Portable Data for Imaging (PDI)  | Portable Media Creator<br>Portable Media Importer<br>Image Display<br>Display | Web Content<br>None<br>None<br>None  |
| <b>Link to Carestream Health IHE information:</b> <a href="http://ihe.carestreamhealth.com">ihe.carestreamhealth.com</a> |   |  |
| <b>Links to Standards Conformance Statements for the Implementation</b>  |   |  |
| <b>DICOM</b>   | <a href="http://dicom.carestreamhealth.com">dicom.carestreamhealth.com</a>    |  |
| <b>Links to general information on IHE</b>   |   |  |
| In North America:<br><a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>  | In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>         | In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a> |