




IHE Integration Statement CARESTREAM Vue PACS v12.0

**Part # 6K9403
2014-10-08**

Integrating the Healthcare Enterprise (IHE) Integration Statement			
Vendor	Product Name	Version	Date
Carestream Health	CARESTREAM Vue PACS: Client Suite, Workflow Manager, IMS, CD Direct Suite, XDS Registry.	CARESTREAM Vue PACS v12.0.0.0	2014- 10-08
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed in this table.			
Note: This document is intended for use by Carestream customers.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None	
Access to Radiology Information (ARI)	Image Manager/Image Archive Image Display Report Repository External Report Repository Access	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display	None	
Cross-enterprise Document Reliable Interchange (XDR)	Document Source	None	
Cross-enterprise Document Sharing for Imaging (XDS.b-l)	Imaging Document Source	Set of DICOM Instances, Multiple documents None	
Cross-enterprise Document Sharing (XDS.b)	Document Repository Document Registry	Multiple document	
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	None	

(PIX)		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Evidence Documents (ED)	Image Manager/Image Archive	None
Simple Image and Numeric Report (SINR)	Report Repository External Report Repository Access	None
Mammography Image (MAMMO)	Image Manager/Image Archive Image Display	None
Key Image Note (KIN)	Image Manager/Image Archive	None
Nuclear Medicine Imaging (NMI)	Image Manager/Image Archive	None
Audit Trail and Node Authentication (ATNA)	Secured Application Image Manager/Image Archive Report Repository Image Display XDS-I Document source	None
Consistent Time (CT)	Time Client	None
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web Content None None None
Link to Carestream Health IHE information: ihe.carestreamhealth.com		
Links to Standards Conformance Statements for the Implementation		
DICOM	dicom.carestreamhealth.com	
Links to general information on IHE		
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j



Carestream Health France
1, rue Galilée
93192 NOISY-LE-GRAND CEDEX
FRANCE



Carestream



Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608
United States

© Carestream Health, Inc. 2014.

Created in the USA.

CARESTREAM is a trademark of Carestream Health, Inc.

Pub No. 6K9403
Rev C