



## Superb Film Duplicates with a Choice of Processing Cycles

Provide referring physicians, insurers, attorneys and others with superb film duplicates that deliver a high level of accurate information. CARESTREAM RA Duplicating Film produces high-quality duplicate radiographic images with virtually one-to-one detail fidelity. Excellent reproduction of the entire sensitometric curve makes this film especially valuable for duplicating chest radiographs.

### Processing Versatility

CARESTREAM RA Duplicating Film can be processed with the same cycles and chemicals used for diagnostic films, to help simplify inventory and processing operations.

Achieve excellent results with your choice of 90- or 60-second processing cycles or their equivalents. Recommended processing chemicals include CARESTREAM RP X-OMAT Developer and Replenisher, CARESTREAM X-OMAT EX II Developer and Replenisher, CARESTREAM X-OMAT LE+ Developer and Replenisher, and CARESTREAM RP X-OMAT LO Fixer and Replenisher.

### Duplicate Radiographs and other Specialty Film Images

This blue-base single-emulsion film is designed for same-size reproduction of radiographs, mammograms and other diagnostic film images.

### Product Information

Film Sizes	CAT No.
35 x 43 cm	837 3805
35 x 35 cm	195 8156
24 x 30 cm	149 4418
18 x 24 cm	114 1159
11 x 14 in	827 5471
10 x 12 in	809 5770
8 x 10 in	811 9117

### Printer Optimization

Duplicating printers must be optimized for use with high-performance CARESTREAM RA Duplicating Film.

This relatively simple procedure can be performed onsite by your authorized Carestream Health representative or distributor/dealer.

*Optimized printers cannot be used with non-Carestream duplicating films.*