

IHE Integration Statement		Date:	February 2008
Vendor	Product Name	Version	
Carestream Health Inc.	<i>Kodak DryView 6800 Laser Imager</i>	1.02	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server	None	
Mammography Image (MAMMO)	Print Server	None	
Internet address for Carestream Health IHE information: http://www.carestreamhealth.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestreamhealth.com/dicom.html		
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe: http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	

End of Document