

IHE Integration Statement		Date:	15 Aug 2003
Vendor	Product Name	Version	
Eastman Kodak Company	<i>Kodak DirectView PACS System 5</i>	5.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Display	None	
Access to Radiology Information (ARI)	Image Display Image Manager/Image Archive	None	
Internet address for Kodak IHE information: www.kodak.com/go/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	www.kodak.com/go/dicom		
Links to general information on IHE			
North America: www.rsna.org/IHE	Europe: www.ihe-europe.org	Japan: www.jira-net.or.jp/ihe-j	