




**IHE Integration Statement
CARESTREAM Vue Cardio PACS
v3.3 Build 7**

**Part # 6M1625
2014-10-08**

IHE Integration Statement			
Vendor	Product Name	Version	Date
Carestream Health	CARESTREAM Vue Cardio PACS	3.3 Build 7	2014-10-04
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed in this table.			
Note: This document is intended for use by Carestream customers.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive	None	
Cardiac Catheterization Workflow (CATH)	Image Manager/Image Archive Performed Procedure Step Manager	Cardiac CATH Intermittently Connected Modality	
Echo Workflow (ECHO)	Image Manager/Image Archive Performed Procedure Step Manager	Echocardiography intermittently Connected Modality	
Simple Image and Numeric Report (SINR)	Report Creator	None	
Consistent Time (CT)	Time Client	SNTP	
Link to Carestream Health IHE information: ihe.carestreamhealth.com			
Links to Standards Conformance Statements for the Implementation			
DICOM	dicom.carestreamhealth.com		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	



Carestream Health France
1, rue Galilée
93192 NOISY-LE-GRAND CEDEX
FRANCE



Carestream



Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608
USA

© Carestream Health, Inc., 2014

Created in the USA.

CARESTREAM is a trademark of Carestream Health.

Pub No. 6M1625
Rev 1.0